

Product datasheet for **RG214464**

IL4R (NM_000418) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	IL4R (NM_000418) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	IL4R
Synonyms:	CD124; IL-4RA; IL4RA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG214464 representing NM_000418
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGGGTGGCTTTGCTCTGGGCTCCTGTCCCTGTGAGCTGCCTGGTCTGCTGCAGGTGGCAAGCTCTG
 GGAACATGAAGGTCTTGCAGGAGCCACCTGCGTCTCCGACTACATGAGCATCTCTACTTGCAGTGGAA
 GATGAATGGTCCCACCAATTGCAGCACCAGCTCCGCCTGTTGTACCAGCTGGTTTTTCTGCTCTCCGAA
 GCCCACACGTGTATCCCTGAGAACACGGAGGCGGGGTGCGTGTGCCACCTGCTCATGGATGACGTGG
 TCAGTGGGATAACTATACTGGACCTGTGGGCTGGGCAGCAGCTGCTGTGGAAGGGCTCCTTCAAGCC
 CAGCGAGCATGTGAAACCCAGGGCCCCAGGAACTGACAGTTCACACCAATGTCTCCGACTCTGCTG
 CTGACCTGGAGCAACCCGTATCCCCTGACAATTACCTGTATAATCATCTCACCTATGCAGTCAACATTT
 GGAGTGAAACGACCCGGCAGATTTGAGATCTATAACGTGACCTACCTAGAACCCTCCCTCCGCATCGC
 AGCCAGCACCCCTGAAGTCTGGGATTTCTACAGGGCACGGGTGAGGGCTGGGCTCAGTGTATAACACC
 ACCTGGAGTGAGTGGAGCCCCAGCACAAGTGGCAACAACCTCTACAGGGAGCCCTTCGAGCAGCACCTCC
 TGCTGGGCGTCAGCGTTTCTGCATTGTATCCTGGCCGTCTGCCTGTTGTGCTATGTACAGCATACCAA
 GATTAAGAAAGAATGGTGGGATCAGATCCCAACCCAGCCCGCAGCCGCCTCGTGGCTATAATAATCCAG
 GATGCTCAGGGGTACAGTGGGAGAAGCGGTCCCGAGGCCAGGAACAGCCAAGTGGCCACACTGGAAGA
 ATTGTCTTACCAAGCTCTTGCCCTGTTTTCTGGAGCACAACATGAAAAGGGATGAAGATCCTCACAAGGC
 TGCCAAAGAGATGCCCTTCCAGGGCTCTGAAAATCAGCATGGTGCCAGTGGAGATCAGCAAGACAGTC
 CTCTGGCCAGAGAGCATCAGCGTGGTGCATGTGTGGAGTTGTTGAGGCCCGGTGGAGTGTGAGGAGG
 AGGAGGAGGTAGAGGAAGAAAAAGGGAGCTTCTGTGCATCGCTGAGAGCAGCAGGGATGACTTCCAGGA
 GGAAGGGAGGGCATTGTGGCCCGCTAACAGAGAGCCTGTTCCCTGGACCTGCTCGGAGAGGAGAATGGG
 GGCTTTTGCAGCAGGACATGGGGGAGTCATGCCTTCTCCACCTTCGGGAAGTACGAGTGTACATGC
 CCTGGGATGAGTCCCAAGTGCAGGGCCCAAGGAGGCACCTCCCTGGGGCAAGGAGCAGCCTCTCCACCT
 GGAGCCAAGTCTCCTGCCAGCCGACCCAGAGTCCAGACAACCTGACTTGCACAGAGACGCCCTCGTC
 ATCGCAGGCAACCCTGCTTACCGCAGCTTACGCAACTCCCTGAGCCAGTACCCTGTCCCAGAGAGCTGG
 GTCCAGACCCACTGCTGGCCAGACACTGGAGGAAGTAGAACCCGAGATGCCCTGTGTCCCCAGCTCTC
 TGAGCCAACCACTGTGCCCAACCTGAGCCAGAAAACCTGGGAGCAGATCCTCCGCCGAAATGTCTCCAG
 CATGGGGCAGCTGCAGCCCCGTCTCGGCCCCACCAGTGGCTATCAGGAGTTTGTACATGCGGTGGAGC
 AGGGTGGCACCCAGGCCAGTGGGTGGTGGGCTTGGGTCCCCAGGAGAGGCTGTTTACAAGGCCTTCTC
 AAGCCTGCTTGCAGCAGTGTGTGTCCCAGAGAAATGTGGGTTTGGGGCTAGCAGTGGGGAAGAGGGG
 TATAAGCCTTCCAAGACCTCATTCTGGCTGCCTGGGGACCCTGCCCCAGTCCCTGTCCCTTGTTC
 CCTTTGGACTGGACAGGGAGCCACCTCGCAGTCCGAGAGCTCACATCTCCAAGCAGCTCCCAGAGCA
 CCTGGGCTGGAGCCGGGGAAAAGGTAGAGGACATGCCAAAGCCCCACTTCCCAGGAGCAGGCCACA
 GACCCCTTGTGGACAGCCTGGGCAGTGGCATTGTCTACTCAGCCCTTACCTGCCACCTGTGCGGCCACC
 TGAAACAGTGTATGGCCAGGAGGATGGTGGCCAGACCCCTGTATGGCCAGTCTTGTGTGGCTGCTG
 CTGTGGAGACAGTCTCGCCCCCTACAACCCCTGAGGGCCCCAGACCCCTCTCCAGGTGGGGTTCCA
 CTGGAGGCCAGTCTGTGTCCGGCTCCCTGGCACCCCTCGGGCATCTCAGAGAAGAGTAAATCCTCATCAT
 CCTTCCATCTGCCCTGGCAATGCTCAGAGCTCAAGCCAGACCCCAAAATCGTGAACCTTGTCTCCGT
 GGGACCCACATACATGAGGGTCTCT

ACCGTACGCGGCCGCTCGAG – GFP Tag – GTTAA

Protein Sequence: >RG214464 representing NM_000418
 Red=Cloning site Green=Tags(s)

```

MGWLCGSLFFPVSCVLLQVASSGNMKVLQEPTCVSDYMSISTCEWKMGPTNCSTELRLLYQLVFLSE
AHTCIPENNGGAGCVCHLLMDDVVSADNYTLDLWAGQQLLWKGSEHVKPRAPGNLTVHTNVSDTLL
LTWSNPYPPDNYLYNHLTYAVNIWSENDPADFRIYNVTYLEPSLRIAASLTKSGISYRARVRAWAQCYN
TWSEWSPSTKWHNSYREPFEQHLLLVSVSCIVILAVCLLCYVSIKIKKEWWDQIPNPARSRLVAIIIQ
DAQSQWEKRSRQEPKCPHWKNCLTKLLPCFLEHNMKREDEPHKAAKEMPFQSGKSAWCPVEISKTV
LWPESISVVRVVELFEAPVECEEEEEVEEEKGFCASPESSRDFQEGREGIVARLTESLFLDLLGEENG
GFCQQDMGESCLLPPSGSTSAHMPWDEFPSAGPKEAPPWGKEQPLHLEPSPPASPTQSPDNLTCTETPLV
IAGNPAYRSFSNSLSQSPCPRELGPDP LLARHLEVEPEMPCVPQLSEPTTVPQPEPETWEQILRRNLQ
HGAAAAPVSAPTSQYQEFVHAVEQGGTQASAVVGLGPPGEAGYKAFSSLLASSAVSPEKCGFGASSGEEG
YKPFQDLIPGCPGDPAPVPVPLFTFGLDREPPRSQSSHLPSSSPEHLGLEPGEKVEDMPKPLPQEAT
DPLVDSLGSIVYSALTCHLCGHLKQCHGQEDGGQTPVMASPCCGCCGDRSSPPTPLRAPDPSPGGVP
LEASLCPASLAPSGISEKSKSSSFHPAPGNAQSSSQTPKIVNFVSVGPTYMRVS
  
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_000418.4](#)

RefSeq Size: 3678 bp

RefSeq ORF: 2478 bp

Locus ID: 3566

UniProt ID: [P24394](#)


Cytogenetics: 16p12.1

Protein Families: Druggable Genome, Secreted Protein

Protein Pathways: Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, Jak-STAT signaling pathway

Gene Summary:

This gene encodes the alpha chain of the interleukin-4 receptor, a type I transmembrane protein that can bind interleukin 4 and interleukin 13 to regulate IgE production. The encoded protein also can bind interleukin 4 to promote differentiation of Th2 cells. A soluble form of the encoded protein can be produced by proteolysis of the membrane-bound protein, and this soluble form can inhibit IL4-mediated cell proliferation and IL5 upregulation by T-cells. Allelic variations in this gene have been associated with atopy, a condition that can manifest itself as allergic rhinitis, sinusitis, asthma, or eczema. Polymorphisms in this gene are also associated with resistance to human immunodeficiency virus type-1 infection. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Apr 2012]

Product images:


Circular map for RG214464