

Product datasheet for **RG214471**

FAM82A1 (RMDN2) (NM_144713) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	FAM82A1
Synonyms:	BLOCK18; FAM82A; FAM82A1; PRO34163; PYST9371; RMD-2; RMD2; RMD4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence: >RG214471 representing NM_144713
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGGAAAGTGCTTAAGCTGCTGCAAAGAGGACCAGAGCTTCCAACGATGTTCTCCAGAAGATCAAGTTA
 GTACAGACGCCAGCATCGTGGAGCCTCTTCTATTTCACAACCAAGTATATCTCTTGGTCATAAAACATC
 TTATTTCCCCTGTAACACATAAAGTCAATGCAGCCAAAGCAAGTAGAAGGTTATTATCTGTATCAAGTCCA
 TCTTTCTCTGAAAGAAGATATTCTTTATTTGTTGGATTTCAAAAAAGAAATGCTTCCCCTTACTGGCAAC
 AAAGTAGAGCTAACTTTGATTCTGAAGAAGACACAGGCTTCACTGATATAAAATCTTCTCAGACCATTG
 TGGTAGCTTTTATTTCCCGCAGAAGACGTTTCTCATCCCGTAAACTAAGTATAGTTTCTTATTACAAGAGT
 GCCATTTTTTTTGTATCCTCAAGCAAGTGGCCGAATGTGTTTAATCTAAATGAAATCGAAATCTTTTCTA
 AAATCTCAAGTAATACTGATGCTAAAAACATATAACCATCTCTGCTCCTGAATATAACACTAAAAATTT
 TAAGAATTTTGAACAAACACTACTTCTCCAGCCTTTGGGAACACTATTGATACAGCCTCCTATCAACAA
 AGCACATCATCCTTCTTTTCTCTTCAAGTGATATTTCTCTCCCGATCAACAAAATGGAATTGCCAATG
 ATATTCACAAAGGGGCCAATTATGTAAGATTTAAAAGATTTCTTTCATCCTCGTCTGAAAGTTACAG
 CACAGATCATTCTCCAATCATGATTCCACAGCATCCCTCTCAAAGTGGAAGTTTCCCATTCCTCCATAAA
 GCTGGATTTTCTTCATCTTATAAAATCTGGTTGCTTTATCCACCTCAAAGCGAATTAAGTTTCAAGGCT
 TGTTTGAAGATGAAGACTTTGCTGTCTTGTCTTCAAGATGAAGACAGATCTTCACTATTGAAATTCCTAA
 AATAAGGTATATTACAGCTAATACTGACACAGAAGAAGAGAGTTTCCAGTCCCTAAGGCATTTAACACA
 CGTGTAGAGGAATTAATTTAGATGCTCTTCTCAGAAGGTAGATCATTACGTATGAGTGAGTCTGGCA
 AGTCGGAGAGTTTGAAGTACTTCTGACACAAAGAAAAGTTTAGAGATGAAATAGAGTTTATGTGGCG
 ATTTGCTCGTGCTTATGGAGACATGTATGAATCTCTACAAACACACAAGAAAAGAAACATTATGCTAAT
 ATTGAAAAAAGTTTAAAGTGAAGAGCTATTAATAGAGCACCATGAATGGACATTGTCATCTGTGGTATG
 CAGTTTTGTGTGGCTATGTATCTGAGTTTGAGGGTTTACAAAACAAAATCAACTATGGGCACCTCTTCAA
 GGAACATCTAGATATAGCAATCAAACTTTTACCAGAGGAACCTTTCTATATTACCTCAAAGGGAGATAC
 TGCTATACTGTCTCAAACTGAGCTGGATTGAGAAAAAATGGCTGCTACTCTGTTTGGAAAAATACCAT
 CTTCAACTGTACAAGAAGCTTTACACAATTTCTTAAGGCTGAAGAAGTATGCCCTGGTTATTCTAATCC
 CAATTACATGTAAGTGTATGTAATGATCTTGAAGGAAAACAGAAATGCTTTGAAGTTCTGTAAT
 TTGGCTTTATTGCTTCTACTGTTACCAAAGAGGAGCTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG214471 representing NM_144713
 Red=Cloning site Green=Tags(s)

MGKCLSCCKEDQSFQRCSPEDQVSTDAQHRGASSISQPSISLGHKTSYSPVTHKVNAAKASRRLLSVSSP
 SFSERRYSLFVGFQKRNASPYWQQSRANFDSEEDTGFTDIKSSSDHCGSFISSRRRFSSRKLISIVSYKS
 AIFFDPPQASGQNVFNLEIEIFSKTSSNTDAKKHITISAPYNTKNFKNFETNTTSPAFGNTIDTASYQQ
 STSSFFSLASDISSPDQQNGIANDIQQRGQLCKDLKDFLHPRPESYSTDHSPIMIPQHPSQSGTFPFLHK
 AGFSSSYKNSGCFIPPQSELTSGLFEDEFVLVFDQEDRSSPIEIPKIRYITANTDTEEQSFVPKAFNT
 RVEELNLDVLLQKVDHLRMSESGKSEFELLRDHKEKFRDEIEFMWRFARAYGDMYELSTNTQEKHYAN
 IGKTLSERAINRAPMNGHCHLWYAVLCGYVSEFEGQLNKINYGHLFKEHLDIAIKLLPEEPFLYYLKGRY
 CYTVSKLSWIEKKMAATLFGKIPSSSTVQEALHNFLKAEELCPGYSNPNYMLAKCYTDLEENQNALKFCN
 LALLLPVTKEEL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_144713

ORF Size: 1719 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_144713.5](#)

RefSeq Size: 1986 bp


RefSeq ORF: 1722 bp

Locus ID: 151393

UniProt ID: [Q96LZ7](#)

Cytogenetics: 2p22.2

Product images:



Circular map for RG214471