

Product datasheet for **RG214666**

ZNHIT3 (NM_004773) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ZNHIT3 (NM_004773) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: ZNHIT3
Synonyms: Hit1; PEHO; TRIP3
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG214666 representing NM_004773
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCGTCGCTCAAATGTAGCACCGTCGCTGCGTGATCTGCTTGGAGAAGCCAAATACCGCTGTCCAG
 CCTGCCGCGTGCCTACTGCTCGGTAGTCTGCTCCGGAAGCACAAGAAGACAGTCAACCCTGAACTCG
 TCCTGTTGAGAAAAAATAAGATCAGCTCTTCTACCAAACCGTAAAGCCTGTGAAAACAAAGATGAT
 GATGACTCTATAGCTGATTTTCTCAATAGTGATGAGGAAGAAGACAGAGTTTCTTGCAGAATTTAAGA
 ATTTAGGGGAATCTGCAACATTAAGAAGCTTATTGCTCAATCCACACCTCAGGCAGTTGATGGTCAACCT
 CGATCAGGGAGAAGACAAAGCAAAGCTCATGAGAGCTTACATGCAAGAGCCTTTGTTTGGAGTTTGCA
 GACTGCTGTTTAGGAATTGTGGACCATCCCAGAATGAGGAGTCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG214666 representing NM_004773
 Red=Cloning site Green=Tags(s)

MASLKCSTVVCVICLEKPKYRCPACRVPYCSVVCFRKHKEQCNPETRPVEKKIRSALPTKTKPVENKDD
 DDSIADFLNSDEEEDRVSLQNLKNGESATLRSLLLPHLRQLMVNLDQGEDKAKLMRAYMQEPLFVEFA
 DCCLGIVEPSQNEES

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_004773

ORF Size: 465 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_004773.4](#)

RefSeq Size: 950 bp

RefSeq ORF: 468 bp


Locus ID: 9326

UniProt ID: [Q15649](#)

Cytogenetics: 17q12

Protein Families: Druggable Genome, Transcription Factors

Product images:



Circular map for RG214666