

Product datasheet for **RG214721**

B3GALT5 (NM_033171) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	B3GALT5 (NM_033171) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	B3GALT5
Synonyms:	3-GalTase 5; B3GalT-V; B3GalTx; B3T5; beta-1; beta-3-Gx-T5; beta3Gal-T5; GLCT5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG214721 representing NM_033171 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTTCCCGAAGATGAGATTGATGTATATTTGCCTTCTGGTTCTGGGGGCTCTTTGTTTGTATTTTA
GCATGTACAGTCTAAATCCTTTCAAAGAACAGTCCTTTGTTTACAAGAAAGACGGAACTCCTTAAGCT
CCCAGATACAGACTGCAGGCAGACACCTCCCTTCCTCGTCCTGCTGGTGACCTCATCCACAAACAGTTG
GCTGAGCGCATGGCCATCCGGCAGACGTGGGGAAAGAGAGGATGGTGAAGGGAAAGCAGCTGAAGACAT
TCTTCCTCCTGGGACCACCAGCAGTGCAGCGGAAACGAAAGAGGTGGACCAGGAGCCAGCGACACGG
GGACATTATCCAGAAGGATTTCCCTAGACGTCTATTACAATCTGACCCTGAAGACCATGATGGGCATAGAA
TGGGTCCATCGCTTTTGTCTCAGGCGGCGTTTGTGATGAAAACAGACTCAGACATGTTCAATGTTG
ACTATCTGACTGAACTGCTTCTGAAGAAAAACAGAAACACCAGGTTTTTCACTGGCTTCTTGAACTCAA
TGAGTTTCCCATCAGGCAGCCATTCAGCAAGTGGTTTGTGAGTAAATCTGAATATCCGTGGGACAGGTAC
CCACCATCTGCTCCGGCACCGGCTACGTGTTTCTGGCGACGTGGCGAGTCAGGTGTACAATGTCTCCA
AGAGCGTCCCACATTAACCTGGAAGACGTGTTTGTGGGGCTCTGCCTCGAAAGGCTGAACATCAGATT
GGAGGAGCTCCACTCCCAGCCGACCTTTTTCCAGGGGGCTTACGCTTCCCGTATGCCCTTTCAGGAGG
ATCGTGGCCTGCCACTTCATCAAGCCTCGGACTCTTTGGACTACTGGCAGGCTCTAGAGAATCCCGGG
GGGAAGATTGTCCGCCTGTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG214721 representing NM_033171
 Red=Cloning site Green=Tags(s)

MAFPKMRMLMYICLLVLGALCLYFSMYSLNPFKEQSFVYKKGDNFLKLPDTCRQTPPFLVLLVTSSHKQL
 AERMAIRQTWGKERMVKGKQLKTFLLGTTSSAAETKEVDQESQRHGDI IQKDFLDVYYNLTKTMMGIE
 WVHRFCPQAAFVMKTDSDMFINVDYLTELLKKNRTRFFTFGLKLNFFPIRQPFKWFVSKSEYPWDRY
 PPFCSGTGYVFSGDVASQYVNVSKSVPIYIKLEDVVFGLCLERLNIRLEELHSQPTFFPGGLRF SVCLFRR
 IVACHFIKPRTLLDYWQALENSRGEDCPPV

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_033171

ORF Size: 930 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_033171.3
RefSeq Size:	2871 bp
RefSeq ORF:	933 bp
Locus ID:	10317
UniProt ID:	Q9Y2C3
Cytogenetics:	21q22.2
Protein Families:	Transmembrane
Protein Pathways:	Glycosphingolipid biosynthesis - globo series, Glycosphingolipid biosynthesis - lacto and neolacto series, Metabolic pathways
Gene Summary:	This gene encodes a member of a family of membrane-bound glycoproteins. The encoded protein may synthesize type 1 Lewis antigens, which are elevated in gastrointestinal and pancreatic cancers. Alternatively spliced transcript variants using multiple alternate promoters have been observed for this gene. [provided by RefSeq, Sep 2017]