

Product datasheet for **RG214741**

ZNF562 (NM_017656) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	ZNF562
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PSI00010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG214741 representing NM_017656
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTCAGCCTTTGATATGTCCCATGATTTCTTTTCTGCTTAACCTCAGAATGGGAAATACAACCTAGAA
CCAAACGGTCATCACTTCAGCAGGGTTTTTTGAAGAATCAAATATTCAGTGGGATACAAATGCAGACAAG
AAGCTACAGTGGATGGAACTCTGTGAGAATTGTGGAGAGGTCTTCAGTGAACAGTTTTGCCTTAAGACA
CACATGAGAGCTCAGAATGGAGGGAACACTTTTGGGGTAATTGTTATGGAAAAGACAGCATCAGTGTGC
ACAAGGAAGCCTCTATTGGACAAGAACTTTCCAACTTAATCCATGTGGAAAAGTCTTCACCTTAACCTCC
AGGCTTGTGTGCATCTTGAAATCTCAATGGAAGACAACCCTACGAATGTAAGGAATGTGGAAAAGGC
TTTAAGTATTTTGAAGCCTTGATAATCACATGGGAATCCACATTGGAGAGAACTCTGTGAATTTTCAGG
AATGTGAGAGAGCCATCACAATTCCTCACACTTAAAGCAGTGTGTAGCAGTTCATACTGGAAAGAAATC
CGAAAAGACTAAGAACTGTGGGAAATCCTTCACTAATTTTCTCAACTTCTGCACATGCAAAAACCTCAT
AAAGGAGAGAAGTCCTTTGAATGTAAGAATGTGGAAGATCCTTTAGAAATTCCTCATCCTTTAATGTTT
ACATTCAAATTCACACTGGAATAAAACACACAAATGTACGGAATGTGGGAAAGCCTTCACTAGATCAAC
TCACCTTACTCAACATGTAAGAACTCACACTGGAATAAAACCTATGAATGTAAGGAATGTGGCCAAGCC
TTCACCTCAGTACACGGGCTTGCTATACATACGAAATCACACTGGAGAGAAACCTATCAGTGTAAAGG
AATGTGGGAGAGCCTTCAATAGATCTTCAACGCTTACTCAACATAGAAGATTACACAGGAGAGAAGCC
TTATGAATGTGTTGAATGTGGGAAGACCTTCATTACTTCTCCCATCGTAGTAACATTTGAAAACCTCAC
AGTGGAGAAAGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG214741 representing NM_017656
Red=Cloning site Green=Tags(s)

MSAFDMSHDFFLCTSEWEIQPRTKRSSLQQGFLKNQIFTGIQMQRYSYSGWKLCEVSEVFSEQFCLKT
HMQAQNNGNTFEGNCYKDSISVHKEASIGQELSKLNPCGVFTLTPGLAVHLEILNGRQPYECKECGKG
FKYFASLDNHMGIHIGELKCEFCERAITTSSHLKQCVAVHTGKKSEKTKNCGKSFTNFSQLSAHAKTH
KGEKSFECKEGRSFRNSSSFNVHIQIHTGIKPHKCTECGKAFTRSTHLTQHVRTHTGIKPYECKECGQA
FTQYTGLAIHNRHTGEKPYQCKEGRFNRSSLTQHRRHTGEKPYECVCEGKTFITSSHRSKHLKTH
SGER

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_017656

ORF Size: 1062 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_017656.2](#)

RefSeq Size: 5541 bp

RefSeq ORF: 1065 bp

Locus ID: 54811

UniProt ID: [Q6V9R5](#)


Cytogenetics: 19p13.2

Domains: zf-C2H2

Protein Families: Transcription Factors

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG214741