

Product datasheet for **RG214762**

ZNF85 (NM_003429) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	ZNF85
Synonyms:	HPF4; HTF1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



ORF Nucleotide Sequence: >RG214762 representing NM_003429
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGGACCATTGACATTTAGGGATGTGGCCATAGAATTCTCTCTGAAGGAGTGGCAATGCCTGGACACTG
 CACAGCGGAATTTATATAGAAATGTGATGTTAGAGAACTACAGAAACCTGGTCTTCCTGGGTATTACTGT
 TTCTAAGCCAGACCTGATCACTTGTCTGGAGCAAGGGAAGAGGCTGGAGTATGAAGAGACATGAGATC
 ATGGTGGCCAAACCCACAGTTATGTGTTCTCATTTGCCCAAGACCTTTGGCCGGAGCAGAAATATAAAG
 ATTCTTTCCAAAAAGTGACACTGAAAAGATATGAAAAATGTAGACATGAAAATTTACCATTAAGAAAAGG
 CTGTGAAAGTATGGATGAGTGTAAAGTGCACAAAGGAGGTTGTAATGGACTTAACCAATGTCTCACAGCT
 ACCCAGAGCAAAATATTCAATGTGATAAATATGTAAAAGTCGCTCATAAATTTCAAATTCAAACAGAC
 ATGAGATAAGACATACTAAAAAGAAACCTTTCAAATGTACAAAATGTGGCAATCATTTGGCATGATTTT
 ATGCCCTAATGAACATAGCAGAAATTCATACTAGAGTAAATTTTACAAAATGTGAAGAATGTGAAAAGCC
 TTTAACTGGTCTCAACCTTACTAAACATAAGAGAATTCATACGGGAGAGAAACCTTACAAATGTGAAG
 AATGTGGTAAAGCCTTTAACCAGTCTCAAACCTTATTAACATAAGAAAATTCATACTGGAGAGAAACC
 CTACAAATGTGAAGAATGTGGCAAACTTTTAACCGATTCTCAACTCTTACTACCCATAAGATAATTCAT
 ACTGGAGAGAAACCTACAAATGTAAAGAATGTGGTAAAGCTTTTAACCGATCTTCAACCTTACTACCC
 ATAGAAAAATTCATACTGGAGAGAAACCTTACAAATGTGAAGAATGTGGCAAAAGCCTTTAAGCAGTCCTC
 AAACCTTACTACACATAAGATAATTCATACTGGAGAGAAACCTACAAATGTAAAAATGTGAAAAGCC
 TTTAACCAGTCTGCACACCTTACCACATGAGGTAATTCATACTGGAGAGAAACCTACAAATGTGAAA
 AATGTGAAAAGCCTTTAATCATTTCTCACACCTTACTACACATAAGATAATTCATACTGGAGAGAAACC
 TTACAAATGTAAAGAATGTGGTAAAGCTTTTAAACACTCTTCAACCTTACTAAACATAAGATAATTCAT
 ACTGGAGAGAAAGCCTTACAAATGTAAAGAATGTGAAAAGCTTTTAACCAATCTCAAACCTTACTGAAC
 ATAAGAAAATTCATACTGGAGAGAAACCTATGAATGTGAAAAATGTGGCAAGCTTTTAACCGATCCTC
 AAATCTTACTAGACATAAGAAAAGTACACAGAGAGAAACCTTACAAATGTGAAGAATGTGGCAAGGT
 TTTAAATGGCCCTCAACCTTACTATCCATAAGATAATTCATACTGGAGAGAAACCTACAAATGTGAAG
 AATGTGGCAAGCTTTTAACCAATCTCAAACCTTACCAACATAAGAAAATTCATACTGGAGAGAAACC
 CTACACATGTGAAGAATGTGGCAAGCCTTTAACCAGTCTCAAACCTTACTAAACATAAGAGAATTCAT
 ACTGGAGAGAAACCTTACAAATGTGAAGAATGTGACAAAGCTTTTAAATGGTCCTCAGTCCTTACTAAAC
 ATAAGATAATTCATACCGAGAGAAAATTACAAATA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG214762 representing NM_003429
 Red=Cloning site Green=Tags(s)

MGPLTFRDVAIEFSLKEWQCLDTAQRNL YRNVML ENYRNL VFLGITVSKPDLITCLEQGKEAWSMKRHEI
 MVAKPTVMCSHFAQDLWPEQNIKDSFQKVTLKRYGKCRHENLPLRKGCESMDECKMHKGGCNGLNQCLTA
 TQSKIFQCDKYVVAHKFSNSNRHEIRHTKKKPFKCTKCGKSGFMISCLTEHSRIHTRVNFYKCEECGKA
 FNWSSTLTCHKRIHTGEKPYKCEECGKAFNQSSNL IKHKKIHTGEKPYKCEECGKTFRNFSTLTTHKIIH
 TGEKPYKCEECGKAFNRSSTLTTHRKIHTGEKPYKCEECGKAFKQSSNL TTHKIIHTGEKPYKCKKCGKA
 FNQSAHLTTHEVIHTGEKPYKCEECGKAFNHFSLTTHKIIHTGEKPYKCEECGKAFKHSSTLTCHKIIH
 TGEKPYKCEECGKAFNQSSKL TEHKKIHTGEKPYEKECGKAFNQSSNL TRHKKSHTEKPYKCEECGKG
 FKWPSTLTIIHKIIHTGEKPYKCEECGKAFNQSSKL TKHKKIHTGEKPYTCEECGKAFNQSSNL TKHKRIH
 TGEKPYKCEECDAFKWSSVL TKHKIIHTGEKLQI

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM_003429

ORF Size: 1785 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_003429.5](#)

RefSeq Size: 2331 bp

RefSeq ORF: 1788 bp

Locus ID: 7639

UniProt ID: [Q03923](#)


Cytogenetics: 19p12

Domains: KRAB, zf-C2H2

Protein Families: Druggable Genome, Transcription Factors

Gene Summary: May be a transcriptional repressor.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG214762