

Product datasheet for **RG214795**

RBIS (NM_001099672) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
 Product Name: RBIS (NM_001099672) Human Tagged ORF Clone
 Tag: TurboGFP
 Symbol: RBIS
 Synonyms: C8orf59
 Mammalian Cell Selection: Neomycin
 Vector: pCMV6-AC-GFP (PS100010)
 E. coli Selection: Ampicillin (100 ug/mL)
 ORF Nucleotide Sequence: >RG214795 representing NM_001099672
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCAAGAACAATTAAGAGGGCCGAAGTCCAGGAATGTATTTACATAGCCAGCCAAAAAACTTTA
 AGGCTAAAAACAAAGCAAAACCAGTTACCACTAATCTTAAGAAGATAAACATTATGAATGAGGAAAAAGT
 TAACAGAGTAAATAAAGCTTTTGTAAATGTACAAAAGGAAGTGCACATTTGCAAAAAGCATTTCATT
 GAACCTCTGCAGAAAGAACTGATTCCTCAGCAGCGTCATGAAAGCAAACCAGTTAATGTTGATGAAGCTA
 CAAGATTAATGGCTCTGTTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG214795 representing NM_001099672
 Red=Cloning site Green=Tags(s)

MAKNKLRGPKSRNVFHIA SQNFKAKNKAKPVTTNLKKINIMNEEKVNRVNKAFVNVQKELAHFAKSISL
 EPLQKELIPQQRHESKPVNVDEATRLMALL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_001099672

ORF Size: 300 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_001099672.2](#)

RefSeq Size: 963 bp

RefSeq ORF: 303 bp

Locus ID: 401466
UniProt ID: [Q8N0T1](#)
Cytogenetics: 8q21.2
Gene Summary: Trans-acting factor in ribosome biogenesis required for efficient 40S and 60S subunit production.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG214795