

Product datasheet for **RG214796**

OR10S1 (NM_001004474) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	OR10S1
Synonyms:	OR11-279
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG214796 representing NM_001004474
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCTGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGACTAGCCGCTCTGTGTGTGAGAAGATGACCATGACAACGGAGAACCCCAACCAGACTGTGGTGAGCC
ACTTCTTCCTGGAGGGTTTGAGGTACACCGCTAAACATTCTAGCCTCTTCTTCCTCTTCTCCTCAT
CTACAGCATCACTGTGGCTGGGAATCTCCTCATCCTCCTAACTGTGGGCTCTGACTCTCACCTCAGCTTA
CCCATGTACCACTTCTGGGGCACCTCTCCTTCTGGATGCCTGTTGTCTACAGTGACAGTGCCCAAGG
TCATGGCAGGCCTGCTGACTCTGGATGGGAAGGTGATCTCCTTTGAGGGCTGTGCCGTACAGCTTTATTG
CTTCCACTTTCTGGCCAGCACTGAGTGCTTCTGTACACAGTCATGGCCTATGACCGCTATCTGGCTATC
TGTC AACCCTGCACTACCCAGTGGCATGAACAGAAGGATGTGTGCAGAAATGGCTGGAATCACCTGGG
CCATAGGTGCCACGCACGCTGCAATCCACACCTCCCTCACCTTCCGCTGCTCTACTGTGGCCTTGCCA
CATTGCCTACTTCTTCTGCGACATACCCCTGTCTAAAGCTCGCCTGTACAGACACCACCAATTAATGAG
CTAGTCATGCTTGCCAGCATTGGCATCGTGGCTGCAGGCTGCCTCATCCTCATCGTTATTTCTACATCT
TCATCGTGGCAGCTGTGTTGCGCATCCGCACAGCCAGGGCCGGCAGCGGGCCTTCTCCCTGCACTGC
CCAGCTCACTGGGTGCTCCTGTACTACGTGCCACCTGTCTGTATCTACCTGCAGCCTCGCTCCAGTGAG
GCAGGAGCTGGGGCCCTGCTGTCTTCTACACAATCGTAACTCCAATGCTCAACCCATTCAATTACACTT
TGCGGAACAAGGAGGTGAAGCATGCTCTGCAAAGGCTTTTGTGCAGCAGCTTCCGAGAGTCTACAGCAGG
CAGCCACCCCA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG214796 representing NM_001004474
Red=Cloning site Green=Tags(s)

MTSRSVCEKMTMTTENPNQTVVSHFFLEGLRYTAKHSSLFFLLFLLIYSITVAGNLLILLTVGSDSHLSL
PMYHFLGHL SFLDACLSTVTPKVMAGLLTLDGKVISFEGCAVQLYCFHFLASTECFLYTMAYDRYLAI
CQPLHYPVAMNRRMCAEMAGITWAI GATHAAIHTSLTFRLLYCGPCHIAFFCDIPPVKLACTDTTINE
LVMLASIGIVAAGCLILIVISYIFIVA AAVLRIRTAQGRQRAFSPCTAQLTGVLLEYVPPVCIYLQPRSSE
AGAGAPAVFYTIIVT PMLNPF IYTLRNKEVKHALQRLLCSSFRETAGSPPPP

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001004474

ORF Size: 993 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001004474.1 , NP_001004474.1
RefSeq Size:	996 bp
RefSeq ORF:	996 bp
Locus ID:	219873
UniProt ID:	Q8NGN2
Cytogenetics:	11q24.1
Protein Families:	Transmembrane
Protein Pathways:	Olfactory transduction
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG214796