

Product datasheet for **RG214861**

SPINT1 (NM_003710) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SPINT1 (NM_003710) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SPINT1
Synonyms:	HAI; HAI1; MANSC2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG214861 representing NM_003710
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCCCTCGAGGACGATGGCCCGGCCCGCTCGCCCGGCCGGCATCCCTGCCGTCGCTTGTGGC
 TTCTGTGCACGCTCGCCCTCCAGGGCACCCAGGCCGGGCCACCGCCCGCCCTGGGCTGCCCGGG
 AGCCGACTGCCTGAACAGCTTTACCGCCGGGTGCCTGGCTTCGTGCTGGACACCAACGCCTCGGTGAG
 AACGGAGTACCTTCTGGAGTCCCCACCGTGCGCCGGGGTGGGACTGCGTGCAGCCTGCTGCACCA
 CCCAGAAGTCAACTTGGCGCTAGTGGAGTGCAGCCGACCGCGGGGAGGACGCCATCGCCGCTGCTT
 CCTCATCAACTGCCTCTACGAGCAGAATTCGTGTGCAAGTTCGCGCCAGGGAGGGCTTCATCAACTAC
 CTCACGAGGGAAGTGTACCGCTCTACCGCAGTGCAGCCAGGGCTTTGGAGGTCTGGGATCCCCA
 AGGCCTGGGACGAGTACTTGAAGGTACAACCCAGGAACCCCTGGTGTGAAGGATGTGAAAAACAC
 AGATTGGCGCCTACTGCGGGGTGACACGGATGTCAGGGTAGAGAGAAAGACCCAAACAGGTGAAACTG
 TGGGGACTCAAGGAAGGCACCTACTGTTCCAGCTGACAGTACTAGCTCAGACCACCCAGAGGACACGG
 CCAACGTACAGTCACTGTGCTGTCCACCAAGCAGACAGAAGACTACTGCCTCGCATCCAACAAGGTGGG
 TCGCTGCCGGGGCTTTTCCACGCTGGTACTATGACCCACGGAGCAGATCTGCAAGAGTTTCGTTTAT
 GGAGGCTGCTTGGGAACAAGAACAACCTACCTCGGAAGAAGAGTGCATTCTAGCCTGTGGGGTGTGC
 AAGGCCCTCCATGAAAAGGCGCCATCCAGTGTGCTCTGGCACCTGTGAGCCACCCAGTTCGCTGCAG
 CAATGGCTGTGCATCGACAGTTTCTGGAGTGTGACGACACCCCAACTGCCCGACGCTCCGACGAG
 GCTGCCTGTGAAAAATACAGAGTGGCTTTGACGAGCTCCAGCGCATCCATTTCCCAAGTACAAAAGGC
 ACTGCGTGGACCTGCCAGACACAGGACTCTGCAAGGAGAGCATCCCGCGTGGTACTACAACCCCTCAG
 CGAACACTGCGCCCGCTTTACCTATGGTGGTGTACGGCAACAAGAACAACCTTTGAGGAAGAGCAGCAG
 TGCTCGAGTCTGTGCGGGCATCTCCAAGAAGGATGTGTTGGCCTGAGCGGGAAATCCCATCCCCA
 GCACAGGCTCTGTGGAGATGGCTGTGCGAGTGTCTGCTCATCTGCATTGGTGGTGGTAGCCATCTT
 GGGTTACTGCTTCTCAAGAACCAGAGAAAGGACTTCCACGGACACCACCACCACCACCACCACCCTC
 GCCAGCTCCACTGTCTCCACTACCGAGGACACGGAGCACCTGGTCTATAACCACACCACGGCCCTC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG214861 representing NM_003710
 Red=Cloning site Green=Tags(s)

MAPARTMARARLAPAGIPAVALWLLCTLGLQGTQAGPPPAPPGLPAGADCLNSFTAGVPGFVLDTNASVS
 NGATFLESPTVRRGWDCVRACCTTQNCNLALVELQPDREGDAIAACFLINCLYEQNFVCKFAPREGFINY
 LTREVYRSYRQLRTQFGGSGIPKAWAGIDLKVPQPEPLVLKDVENTDWRLLRGDTRVVERKDPNQVEL
 WGLKEGTYLFLQTVTSSDHPEDTANVTVTVLSTKQTEDYCLASNKVGRCRGSFPRWYDPTQICKSFVY
 GGCLGNKNNYLREEECILACRGVQGPSMERRHPVCSGTQPTQFRCSNGCCIDSFLECDTPNCPDASDE
 AAKEYTSGFDELQRIHFPSDKGHCVDLPDTGLCKESIPRWYYPFSEHCFRFTYGGCYGNKNNFEEQQ
 CLESCRGISKKDVFGLRREIPIPSTGSVEMAVAVFLVICIVVVVAAILGYCFFKNQRKDFHGHHPPTP
 ASSTVSTTEDTEHLVYNHTTRPL

TRTRPLE – GFP Tag – V

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:


ACCN: NM_003710

ORF Size: 1539 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_003710.4](#)

RefSeq Size: 2436 bp

RefSeq ORF: 1542 bp

Locus ID: 6692

UniProt ID: [O43278](#)


Cytogenetics: 15q15.1

Domains: Idl_recept_a, KU

Protein Families: Druggable Genome, Secreted Protein, Transmembrane

Gene Summary: The protein encoded by this gene is a member of the Kunitz family of serine protease inhibitors. The protein is a potent inhibitor specific for HGF activator and is thought to be involved in the regulation of the proteolytic activation of HGF in injured tissues. Alternative splicing results in multiple variants encoding different isoforms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG214861