

Product datasheet for **RG214934**

UBE2E1 (NM_182666) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: UBE2E1 (NM_182666) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: UBE2E1
Synonyms: UBCH6
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG214934 representing NM_182666
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTCGGATGACGATTTCGAGGGCCAGCACCAGCTCCTCCTCATCTTCGTCTTCCAACCAGCAAACCGAGA
 AAGAAACAACACCCCAAGAAGAAGGAGAGTAAAGTCAGCATGAGCAAAAACCTCAAACCTCTCCAC
 CAGCGCAAGAGTGTGGTCCCAAAGGCGATAACATCTATGAATGGAGATCAACCATTCTAGGGCTCCA
 GGATCCGTGTATGAGGGTGTATTCTTCTCGATATCACTTTACACCAGAATATCCCTTCAAGCCTC
 CAAAGTTACATTTCCGACAAGAATCTATCATTGTAATATTAACAGTCAAGGTGTTATTTGCTTGGACAT
 ATTGAAAGATAATTGGAGTCCAGCACTAACATTTCTAAAGTCCCTTTCTATCTGCTCACTTCTTACA
 GACTGTAATCCTGCCGCCCTTGGTGGGAAGTATTGCCACTCAGTATATGACCAACAGAGCAGAACATG
 ACAGAATGGCCAGACAGTGGACCAAGAGATACGCTACA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG214934 representing NM_182666
 Red=Cloning site Green=Tags(s)

MSDDDSRASTSSSSSSSNQQTEKETNTPKKKESKVSMSKNSKLLSTSAKSAGPKGDNIYWRSTILGPP
 GSVYEGGVFFLDITFTPEYFPKPPKVTFRTRIIYHCNINSQGVICLDILKDNWSPALTISKVLLSICSLLT
 DCNPADPLVGSIAQYMTNRAEHRMARQWTKRYAT


TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_182666

ORF Size: 528 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_182666.3
RefSeq Size:	1428 bp
RefSeq ORF:	531 bp
Locus ID:	7324
UniProt ID:	P51965
Cytogenetics:	3p24.2
Protein Pathways:	Ubiquitin mediated proteolysis
Gene Summary:	The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. Three alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jan 2011]