

## Product datasheet for **RG214979**

### FGFR1 (NM\_023108) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	FGFR1 (NM_023108) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	FGFR1
Synonyms:	BFGFR; CD331; CEK; FGFBR; FLG; FLJ99988; FLT2; HBGFR; KAL2; N-SAM; OGD
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG214979 representing NM_023108 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGTGGAGCTGGAAGTGCCTCCTCTTCTGGGCTGTGCTGGTCACAGCCACACTCTGCACCGCTAGGCCGT  
CCCCGACCTTGCCTGAACAAGATGCTCTCCCCTCCTCGGAGGATGATGATGATGATGACTCCTCTTC  
AGAGGAGAAAGAAACAGATAACACCAAACCAACCCCGTAGCTCCATATTGGACATCCCCAGAAAAGATG  
GAAAAGAAATTGCATGCAGTGCCGGCTGCCAAGACAGTGAAGTTCAAATGCCCTTCCAGTGGACCCCAA  
ACCCACACTGCGCTGGTTGAAAAATGGCAAGAATTCAAACCTGACCACAGAATTGGAGGCTACAAGGT  
CGTTATGCCACCTGGAGCATCATAATGGACTCTGTGGTGCCTCTGACAAGGGCAACTACACCTGCATT  
GTGGAGAATGAGTACGGCAGCATCAACCACACATACCAGCTGGATGTCGTGGAGCGGTCCCTCACCGGC  
CCATCCTGCAAGCAGGTTGCCCGCAACAAAACAGTGGCCCTGGGTAGCAACGTGGAGTTCATGTGTAA  
GGTGTACAGTGACCCGCAGCCGCACATCCAGTGGCTAAAGCACATCGAGGTGAATGGGAGCAAGATTGGC  
CCAGACAACCTGCCTTATGTCCAGATCTGAAGGTAATCATGGCACCAGTCTTCGTGGGCCAGTCTACTG  
GGAAGGAGACCACTGTCTCGGGGCTCAAGTTCCTGTGGGCAGGCTCAGTTGCCCCGAATGGGATCATT  
CCTCACGCTTCAGGCACACACTCCATCTCAGTAGGGATCTAGCCACATCCCCAGGACTAGTAACAGA  
GGTCACAAAGTGGAGGTGAGCTGGGAACAGAGGGCTGCAGGGATGGGTGGTCTGGTCTG

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG214979 representing NM\_023108  
Red=Cloning site Green=Tags(s)

MWSWKCLLFWAVLVATLCTARPSPTLPEQDALPSEDDDDDDSSSEKETDNTKPNPVAPYWTSPEKM  
 EKKLHAVPAAKTVKFKCPSGTPNPTLRWLKNGKEFKPDHRIGGYKVRYATWSIIMDSVVP SDKGN YTCI  
 VENEYGSINHTYQLDVVERSPHRPILQAGLPANKTVALGSNVEFMCKVYSDPQPHIQWLKHIEVNGSKIG  
 PDNLPYVQILKVI MAPV FVGQSTGKETTVSQAQVPVGR LSCPRMGSFLTLQAHTLHL SRDLATSPRTSNR  
 GHKVEVSWEQRAAGMGGAGL

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



**Plasmid Map:**



**ACCN:** NM\_023108

**ORF Size:** 900 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_023108.2</a> , <a href="#">NP_075596.1</a>
<b>RefSeq Size:</b>	2583 bp
<b>RefSeq ORF:</b>	902 bp
<b>Locus ID:</b>	2260
<b>Cytogenetics:</b>	8p11.23
<b>Domains:</b>	ig
<b>Protein Families:</b>	Druggable Genome, Protein Kinase, Transmembrane
<b>Protein Pathways:</b>	Adherens junction, MAPK signaling pathway, Melanoma, Pathways in cancer, Prostate cancer, Regulation of actin cytoskeleton

**Gene Summary:**

The protein encoded by this gene is a member of the fibroblast growth factor receptor (FGFR) family, where amino acid sequence is highly conserved between members and throughout evolution. FGFR family members differ from one another in their ligand affinities and tissue distribution. A full-length representative protein consists of an extracellular region, composed of three immunoglobulin-like domains, a single hydrophobic membrane-spanning segment and a cytoplasmic tyrosine kinase domain. The extracellular portion of the protein interacts with fibroblast growth factors, setting in motion a cascade of downstream signals, ultimately influencing mitogenesis and differentiation. This particular family member binds both acidic and basic fibroblast growth factors and is involved in limb induction. Mutations in this gene have been associated with Pfeiffer syndrome, Jackson-Weiss syndrome, Antley-Bixler syndrome, osteoglophonic dysplasia, and autosomal dominant Kallmann syndrome 2. Chromosomal aberrations involving this gene are associated with stem cell myeloproliferative disorder and stem cell leukemia lymphoma syndrome. Alternatively spliced variants which encode different protein isoforms have been described; however, not all variants have been fully characterized. [provided by RefSeq, Jul 2008]