

Product datasheet for **RG215112**

GSG1 (NM_001080555) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: GSG1 (NM_001080555) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: GSG1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG215112 representing NM_001080555
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGATCGCC

ATGAGCGATCCCTCTCAACTGACTCAAATGTTTGCCTCACCCAGGAGATGGAGCTCTCGAAGGCCTTCT
 CTGGCCAGCGGACACTCCTATCTGCCATCCTCAGCATGCATCACTCAGCTTCTCCACAACATCCCTGCT
 CAGCAACTACTGGTTTGTGGGCACACAGAAGGTGCCAAGCCCTGTGCGAGAAAGGTCTGGCAGCCAAG
 TGCTTTGACATGCCAGTGTCCCTGGATGGAGATACCAACACATCCACCCAGGAGGTGTTACAATAACAAC
 GGGAGACTGGGGATGACCGTTCTCCTCCGGAGCTCCGGAGTGGCATGTGGCTATCCTGTGAGGAAAC
 TGTGGAAGAACCAGCACTGCTCCATCCCCAGTCCTGGAACAATTTAGAGCCCTTCGGTCCAGTGGTACA
 GCGGCAGCAAAAGGGGAGAGGTGCCGAAGTTTCATTGAACTTACACCACCAGCCAAGAGAGAAATCCTAT
 GGTTATCCCTGGGAACGCAGATCACCTACATCGGACTTCAATTCATCAGCTTCTCCTGCTACTAACAGA
 CTTGCTACTCACTGGGAACCCTGCCTGTGGGCTCAAACCTGAGCGCCTTTGCTGCTGTTTCTCTGCTCCTG
 TCAGGTCTCCTGGGGATGGTGGCCACATGATGTATTACAAGTCTTCCAAGCGACTGTCAAACCTGGGTC
 CAGAAGACTGGAGACCACATGTTTGAATTATGGCTGGGCCTTCTACATGGCCTGGCTCTCCTTACCTG
 CTGCATGGCGTCGGCTGTACCACCTTCAACACGTACACCAGGATGGTGTGGAGTTCAAGTCAAGCAT
 AGTAAGAGCTTCAAGGAAAACCCGAACCTACCACATCACCATCAGTGTTCCTCGGGGCTGTCAA
 GTGCAGCCCCACCGTGGGTCTTTGACCAGCTACCACCAGTATCATAATCAGCCATCCACTCTGTCTC
 TGAGGGAGTGCAGTTCTACTCCGAGCTGCGGAACAAGGGATTTCAAGAGGGGCCAGCCAGGAGCTGAAA
 GAAGCAGTTAGGTCATCTGTAGAGGAAGAGCAGTGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG215112 representing NM_001080555
Red=Cloning site Green=Tags(s)


MSDPSQLTQNVCLTQEMELSKAFSGQRTLLSAILSMLSLSFSTTSLLSNYWVFGTQKVPKPLCEKGLAAK
 CFDMPVSLDGDNTSTQEVVQYNWETGDDRF SFRSFRSGMWLSCEETVEEPALLHPQSWKQFRALRSSGT
 AAKGERCRSFIELTPPAKREILWLSLGTQITYIGLQFISFLLLLDLLLLTGNPACGLKLSAFAAVSSVL
 SGLLGMVAHMMYSQVFQATVNLGPEDWRPHVWNYGWAFYMAWLSFTCCMASAVTTFNTYTRMVLEFKCKH
 SKSFKENPNCLPHHHQCFPRRLSSAAPTVGPLTSYHQYHNQPIHSVSEGVDYSELRNKGFQRGASQELK
 EAVRSSVEEEQC

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_001080555

ORF Size:	1086 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001080555.4
RefSeq Size:	2055 bp
RefSeq ORF:	1089 bp
Locus ID:	83445
Cytogenetics:	12p13.1
Protein Families:	Transmembrane
Gene Summary:	May cause the redistribution of PAPOLB from the cytosol to the endoplasmic reticulum. [UniProtKB/Swiss-Prot Function]