

Product datasheet for **RG215206**

IMPDH1 (NM_000883) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	IMPDH1 (NM_000883) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	IMPDH1
Synonyms:	IMPD; IMPD1; IMPDH-I; LCA11; RP10; sWSS2608
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG215206 representing NM_000883
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGAGGGGCCACTCACTCCACCACCGCTGCAGGGAGGCGGAGCCGCGTGTCCGGAGCCCGGAGCC
 GGCAACACCCGGACACGAGACGGCGCGCAGCGGTACAGCGCCGACTGCTGCAGGCGCGCTACGAGCC
 CGAGAGCCCTAGATTGGACCTCGCTACACACCCGACGACACCCCGTTTCAGAACTATCTTCAGTGGTCTTA
 CTGGCAGGTGTTGGTGTCCAGATGGATCGCCTTCGAGGGCTAGCATGGCGGACTACCTGATCAGCGGCG
 GCACCGGTACTGTGCCGAGGATGGGCTCACCGCGCAGCAGCTCTTCGCCAGCGCCGACGGCCTCACCTA
 CAACGACTTCTGATTCTCCAGGATTCATAGACTTCATAGCTGATGAGGTGGACCTGACCTCAGCCCTG
 ACCCGGAAGATCAGCTGAAGACGCCACTGATCTCTCCCCCATGGACTGTGACAGAGGCTGACATGG
 CCATTGCCATGGCTCTGATGGGAGGATTGGTTTCATTACCACAACGCACCCAGAGTCCAGGCCAA
 CGAGGTGCGGAAGTCAAGAAGTTTGAACAGGGCTTCATCAGGACCCTGTGGTGTGAGCCCTCGCAC
 ACTGTGGGCGATGTGCTGGAGGCCAAGATGCGGCATGGCTTCTCTGGCATCCCCATCACTGAGACGGGCA
 CCATGGGCGCAAGCTGGTGGGCATCGTCACCTCCCGAGACATCGACTTCTTGCTGAGAAGGACCACAC
 CACCCTCCTCAGTGAGGTGATGACGCCAAGGATTGAAGTGGTGGTGGCTCCAGCAGGTGTGACGTTGAAA
 GAGGCAATGAGATCCTGCAGCGTAGCAAGAAAGGGAAGTGCCTATCGTCAATGATTGCGATGAGCTGG
 TGGCCATCATCGCCCGCACCGACTGAAGAAGAACCGAGACTACCTCTGGCCTCCAAGGATTCCAGAA
 GCAGCTGCTCTGTGGGCGAGCTGTGGGCACCCGTGAGGATGACAAATACCGTCTGGACCTGCTCACCAG
 GCGGGCGTCGACGTCATAGTCTTGACTCGTCCAAGGGAATTCGGTGTATCAGATCGCCATGGTGCATT
 ACATCAAAACAGAAGTACCCACCTCCAGGTGATTGGGGGAACGTGGTGACAGCAGCCAGGCCAAGAA
 CCTGATTGATGCTGGTGTGGACGGCTGCGCGTGGGCATGGGCTGCGGCTCCATCTGCATACCCAGGAA
 GTGATGGCCTGTGGTGGCCCGGACTGCTGTGTACAAGGTGGCTGAGTATGCCCGGCGCTTTGGTG
 TGCCCATCATAGCCGATGGCGGCATCCAGACCGTGGGACACGTGGTCAAGGCCCTGGCCCTTGAGCCTC
 CACAGTGTGATGGGCTCCCTGCTGGCCGCCACTACGGAGGCCCTGGCGAGTACTTCTTCTCAGACGGG
 GTGCGGCTCAAGAAGTACCGGGCATGGGCTCACTGGATGCCATGGAGAAGAGCAGCAGCAGCCAGAAAC
 GATACTTCAGCGAGGGGATAAAGTGAAGATCGCRCAGGTGTCTCGGGTCCATCCAGGACAAAGGATC
 CATTGAGAAGTTCGTGCCCTACCTCATAGCAGGCATCCAACACGGCTGCCAGGATATCGGGCCCGCAGC
 CTGTCTGTCTTCGGTCCATGATGTACTCAGGAGAGCTCAAGTTTGAGAAGCGGACCATGTGCGCCAGA
 TTGAGGGTGGTGTCCATGGCCTGCACTTACGAAAAGCGGCTGTAC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG215206 representing NM_000883
 Red=Cloning site Green=Tags(s)

MEGPLTPPPLQGGGAAVPEPGARQHPGHETAARYSARLLQAGYEPESPRLDLATHPTTNRSELSSVVL
 LAGVGVQMDRLRRASMADYLISSGTYVPEDGLTAQQLFASADGLTYNDFLILPGFIDFIADVDLTSAL
 TRKITLKTPLISSPMDTVTEADMAIAMALMGGIGFIHNNCTPEFQANEVRKVKKFEQGFITDPVVLSPSH
 TVGDVLEAKMRHGFSGIPITETGTMGSKLVGIVTSRDIDFLAEKDHTLLSEVMTPRIELVVAPAGVTLK
 EANEILQRSKKGLPIVNDCELVAIARTDLKKNRDYPLASKDSQKQLLCGAAVGTREDDKYRLDLLTQ
 AGVDVIVLDSSQNSVYQIAMVHYIKQKYPHLQVIGGNVVTAAQAKNLIDAGVDGLRVGMGCGSICITQE
 VMACGRPQGTAVYKVAEYARRFVPIIADGGIQTVGHVVKALALGASTVMMGSLLAATTEAPGEYFFSDG
 VRLKKYRGMGSLDAMEKSSSSQKRYFSEGDKVKIXQGVSGSIQDKGSIQKFPYPLIAGIQHGCDIGARS
 LSVLRSMYSSELKFEKRTMSAQIEGGVHGLHSYEKRLY

TRTRPLE – GFP Tag – V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_000883

ORF Size: 1797 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_000883.2](#), [NP_000874.2](#)

RefSeq Size: 2880 bp

RefSeq ORF: 1800 bp


Locus ID: 3614

UniProt ID: [P20839](#)

Cytogenetics: 7q32.1

Domains:	CBS, IMPDH
Protein Families:	Druggable Genome
Protein Pathways:	Drug metabolism - other enzymes, Metabolic pathways, Purine metabolism
Gene Summary:	The protein encoded by this gene acts as a homotetramer to regulate cell growth. The encoded protein is an enzyme that catalyzes the synthesis of xanthine monophosphate (XMP) from inosine-5'-monophosphate (IMP). This is the rate-limiting step in the de novo synthesis of guanine nucleotides. Defects in this gene are a cause of retinitis pigmentosa type 10 (RP10). Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2008]

Product images:



Circular map for RG215206