

Product datasheet for **RG215224**

SIRPG (NM_001039508) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SIRPG (NM_001039508) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	SIRPG
Synonyms:	bA77C3.1; CD172g; SIRP-B2; SIRPB2; SIRPgamma
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG215224 representing NM_001039508 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCTGTCCCAGCCTCCTGGCCCCATCCTCCTGGTCCTTTCCTGCTTCTGACTCTACTGCTGGGACTTA
CAGAAGTGGCAGGTGAGGAGGAGCTACAGATGATTCAGCCTGAGAAGCTCCTGTTGGTCACAGTTGGAAA
GACAGCCACTCTGACTGCACTGTGACCTCCCTGCTTCCCGTGGGACCCGTCCTGTGGTTCAGAGGAGTT
GGACCAGGCCGGGAATTAATCTACAATCAAAAAGAAGGCCACTTCCCAGGGTAACAACAGTTTCAGACC
TCACAAAGAGAAACAACATGGACTTTTCCATCCGCATCAGTAGCATCACCCAGCAGATGTCGGCACATA
CTACTGTGTGAAGTTTCGAAAAGGGAGCCCTGAGAACGTGGAGTTAAGTCTGGACCAGGCACTGAGATG
GCTTTGGGTGCCAAACCCTCTGCCCCGTGGTATTGGGCCCTGCGGCGAGGACCACACCTGAGCATAACAG
TGAGTTTCACCTGTGAGTCCCATGGCTTCTCTCCAGAGACATCACCTGAAATGGTTCAAAAATGGGAA
TGAGCTCTCAGACTTCCAGACCAACGTGGACCCACAGGACAGAGTGTGGCTACAGCATCCGCAGCACA
GCCAGGGTGGTACTGGACCCCTGGGACGTTGCTCTCAGGTCATCTGCGAGGTGGCCCATGTACCTTGC
AGGGGGACCCTCTTCGTGGGACTGCCAATTGTCTGAGGCCATCCGAGGCCCGGCATCATCCCTTACTGC
GCTGCTCTCATAGCTGTCTCCTGGGCCCATCTACGTCCCCTGGAAGCAGAAGACC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG215224 representing NM_001039508
 Red=Cloning site Green=Tags(s)

MPVPASWPHPPGPFLLLTLGLTEVAGEEELQMIQPEKLLLVTVGKTATLHCTVTSLLPVGPVLWFRGV
 GPGRELIYNQKEGHFPRVTTVSDLTKRNNMDFSIRISSITPADVGTYYCVKFRKGSPEVFKSGPGTEM
 ALGAKPSAPVVLGPAARTTPEHTVSFTCESHGFSRPRDITLKWFKNGNELSDFQTNVDPTGQSVAYSIRST
 ARVVLDPWDVRSQVICEVAHVTLQGDPLRGTANLSEAIRGPASSLTALLLIAVLLGPIYVPWKQKT

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_001039508

ORF Size: 828 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001039508.1 , NP_001034597.1
RefSeq Size:	1387 bp
RefSeq ORF:	831 bp
Locus ID:	55423
UniProt ID:	Q9P1W8
Cytogenetics:	20p13
Protein Families:	Druggable Genome, Transmembrane
Gene Summary:	The protein encoded by this gene is a member of the signal-regulatory protein (SIRP) family, and also belongs to the immunoglobulin superfamily. SIRP family members are receptor-type transmembrane glycoproteins known to be involved in the negative regulation of receptor tyrosine kinase-coupled signaling processes. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]