

Product datasheet for **RG215254**

KRTAP10-9 (NM_198690) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KRTAP10-9 (NM_198690) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	KRTAP10-9
Synonyms:	KAP10.9; KAP18.9; KRTAP18-9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG215254 ORF sequence, codon optimized. Due to the complexity of NM_198690, the ORF clone is codon optimized for mammalian Expression. The nucleotide sequence differs from the reference sequence, yet the amino acid sequence remains identical.

Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCGCTTCAACAATGTCCATCCGCTCATCTGCTTACAGCGATTATGGCAGGTGGATGATTGTCCTG
AGAGCTGTTGCGAGCCCCGTGCTGCGCTACCTCCTGCTGTGCCCGGCTCCCTGTCTCACCTGGTCTG
TACACCACTCTCTAGAGTGTCTAGTCCCTGCTGTGAGGTGACGTGCGAACCAAGCCCATGCCAGTCAGGT
TGTACCACTCTCTGCACTCCTTCTTGTGCCAGCAATCCAGCTGCCAGCCCGCTACTGTACCTCTTCCC
CCTGCCAGCAAGCCTGCTGCGTTCCAGTTTGTGTAAACCTGTGTGCTGTGTCCCTGTGTGTTGTGGAGC
GAGCTCATGTTGTCAACAATCCTCTTACCAACCTGCTTGTGTGCTGCTCCAGTAGTTGTGACCCCGCATGC
TGTGTGCCAGTCTGCTGCAAGCCGGTATGTTGCGCGCTACATGCTCCGAAGATAGTTACAGCTGTTGCC
AGCAGTCCAGTTGTCAACCCGCTTGTGCAACAGCTCTCCATGTCAACAGTCTTACTGTGTGCCGGTGTG
CTGTAAGCCAGTGTGTTGTAAACCATCTGTTGCGTACCGTTTGTAGCGCGCGAGCTCACTGTGTTGC
CAACAGTCTGGATGCCAGCCTGCATGCTGCACTACCAGTTGTTGTGCGCCATCCAGCTCCGTGTCCCTGC
TTTGTGCGCCAGTTGTGCAACCCGCTGTTGCGTCCCCGTGCTCTTGTGCGCACCCACTAGCTCTAG
GCAACCCAGCTACTGCCGCGAGGCGAGCTGCGTTAGTCTGCTCTGTGCGCCGTTTGTCTCGACCAAGCC
TGTTACAGTTTTAGTAGCGGGCAGAAGTCATCTTGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG215254 representing NM_198690
 Red=Cloning site Green=Tags(s)

MAASTMSIRSSAYSDSWQVDDCPESCCEPPCCATSCCAPAPCLTLVCTPVSrvSSPCCQVTCEPSPCQSG
 CTSSCTPSCCQQSSCQPAYCTSSPCQQACCVPVCKPVCCVPVCCGASSCCQSSYPACASSSCQPAC
 CVPVCCKPVCCAPTCSSEDSYSCCQSSCQPACCTSSPCQQSYCVPVCKPVCCKPICCVPCSGASSLCC
 QQSGCQPACCTTSCCRPSSSVLLCRPVC RPACCVPVSSCCAPTSSRQPSYCRQASCVSLLCRPVCSRPA
 CYFSFGQKSSC

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_198690

ORF Size: 876 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_198690.2, NP_941963.2

RefSeq Size: 1276 bp


RefSeq ORF: 879 bp

Locus ID: 386676

UniProt ID: P60411

Cytogenetics: 21q22.3

Product images:



Circular map for RG215254