

Product datasheet for **RG215338**

FOXP4 (NM_001012427) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FOXP4 (NM_001012427) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	FOXP4
Synonyms:	hFKHLA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG215338 representing NM_001012427
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGATGGTGAATCTGCCTCGGAGACAATCAGGTCGGCTCCATCTGGTCAGAATGGCGTGGGCAGCCTCT
 CTGGCAAGCCGATGGCAGCAGCGCGGGCCACAGGGACAACGCAAGTGGCACGGGCAGGGAAGTGAC
 CACGGGTGCAGACAGCAATGGTGAGATGAGTCCCGCAGAGCTGCTGCACTTCCAGCAGCAACAGGCTCTC
 CAAGTGGCCCGCAGTTCTGCTGCAGCAGGCCTCAGGCCTGAGCTCCCCAGGGAACAATGACAGCAAAC
 AGTCTGCCTCTGCTGTGACAGGTGCCTGTGTCGGTGGCCATGATGTCGCCGAGATGCTTACCCCGCAACA
 GATGCAGCAGATCCTGTCGCCCCCGCAGCTGCAGGCCTTGTCCAGCAGCAGCAAGCCCTCATGCTCCAG
 CAGGAGTACTACAAGAAGCAGCAGGAGCAGCTCCACCTGCAGCTCCTCACCCAGCAGCAGGCTGGGAAAC
 CGCAGCCCAAAGAGGCACTGGGGAACAAGCAGCTGGCCTTCCAGCAGCAGCTCCTGCAAAATGCAACAGTT
 GCAGCAGCAGCACCTGCTCAACCTGCAGAGGCAGGGGCTGGTCAGCCTGCAGCCCAACCAAGCCTCGGGG
 CCCTCCAGACCTTCCGCAAGCAGCTGTTTGCCCAACAGACCTGCCCCAGCTGTGGAAGGGCGAGGGTG
 CCCCCGGGCAGCCTGCCGAGGACAGCGTCAAGCAGGAGGGGCTGGACCTCACTGGCACGGCCGCCACCGC
 TACCTCGTTTGCCGCTCCCCCAAGGTCTCACCCCTCTCCCACCATACCTGCCCAACGGACAGCCT
 ACTGTGCTCACATCTCGGAGAGACAGCTTCCACGAGGAGACCCCGGCTCCCACCCCTGTACGGAC
 ACGGAGAGTGCAAGTGGCCAGGCTGTGAGACCTGTGTGAAGACCTGGGCCAGTTTATCAAACACCTCAA
 CACAGAGCACGCCCTGGATGACCGGAGTACAGCCAGTCCGGGTACAGATGCAGGTGGTGCAGCAGCTG
 GAGATCCAGCTCGCAAGGAGAGCGAGCGGCTGCAGGCCATGATGGCCACCTGCACATGCGGCCCTCGG
 AGCCCAAGCCCTTCAGCCAGCCACTGAACCCGGTCCCCGGCTCCTCCTCATTCTCCAAGGTGACCCGTCTC
 TGACAGACTCATTCCAGATGGTCTCGTGCACCCCGACCTCGGCCGAGCCCTGTACCCCTCTA
 CGGCCCTTGGCCTGGGCTCTGCCTCCCTGCATGGTGGGGGCCAGCCGTCGGAGAAGCAGTGACAAGT
 TCTGCTCCCCATCTCCTCAGAGCTGGCCAGAATCATGAGTTCTACAAGAACGCCGACGTCCGGCCCC
 CTTACCTACGCCTCCCTCATCCGCCAGGCCATCCTGAAACCCCTGACAGGCAGCTGACCTGAATGAG
 ATCTATAACTGGTTCACCAGGATGTTGCTATTTCCGCAGAAACACTGCCACCTGGAAGAACGCCGTGC
 GCCAACCTCAGCCTGCACAAGTCTGTCGCGTGGAGAACGTCAGGGTGGCGTGTGGACTGTGGA
 CGAGCGGGAGTATCAGAAGCGGAGACCCCAAAGATGACAGGGAGCCCCACCTGGTGAAGAACATGATC
 TCTGGCCTCAGCTATGGAGCACTTAATGCCAGCTACCAGGCCCTGGCCGAGAGCAGCTTCCCCCTCC
 TCAACAGCCCTGGCATGCTGAACCTGGTCCGCCAGCAGCCTGCTGCCCTCAGCCACGATGACGTGGG
 TGCCCCGTGGAGCCGCTGCCAGCAACGGCAGCAGCAGCCCTCCTCGCCTCTCCCCGCCCCAGTACAGC
 CACCAGGTGCAGGTGAAGGAGGAGCCAGCAGAGGCAGAGGAAGACAGGCAGCCCGGGCCTCCCCTGGGCG
 CCCCTAACCCAGCGCTCGGGGCTCCGGAAGACAGGGACCTGGAGGAGGAGCTGCCGGGAGAAGAACT
 GTCC

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG215338 representing NM_001012427
 Red=Cloning site Green=Tags(s)


```

MMVESASETIRSAPSGQNGVGSLSGQADGSSGGATGTTASGTGREVTTGADSNNGEMSPAELLHFQQQAL
QVARQFLLQQASGLSSPGNNDKQASASAVQVPVSVAMMSPQMLTPQQMQQILSPPQLQALLQQQQALMLQ
QEYKKQEQELHLQLLTQQQAGKPQKEALGNKQLAFQQQLLQMQQLQQQHLLNLQRQLVSLQPNQASG
PLQTLPQAAVCPDLPQLWKGEAGPQPAEDSVKQEGLDLTGTAATATSF AAPPKVSPPLSHHTLPNGQP
TVLTSRRDSSHEETPGSHPLYGHGECKWPGCETLCELDLQGIKHLNTEHALDDRSTAQCRVQMQRVQQL
EIQLAKESESLQAMMAHLHMRPSEPFPKFSQPLNVPVPGSSSF SKVTVSAADSFDPGLVHPPTSAAPVTP
RPPGLGSASLHGGPARRRSSDKFCSPISSELAQNHEFYKNADVRPFTYASLIRQAILETDRQLTLNE
IYNWFTRMFAYFRRNTATWKNVVRHNL SLHKCFVRVENVKAVVTVDEREYQRRPPKMTGSPTLVKMI
SGLSYGALNASYQAALAESSFPLLNSPGMLNPGSASSLLPLSHDDVGPVPELPSNGSSSPRLSPPQYS
HQVQVKEEPAEAEEDRQPGPPLGAPNPSASGPPEDRDLEELPGEELS
  
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001012427

ORF Size: 2034 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001012427.2](#)

RefSeq Size: 5959 bp

RefSeq ORF: 2037 bp

Locus ID: 116113


UniProt ID: [Q8IVH2](#)

Cytogenetics: 6p21.1

Protein Families: Transcription Factors

Gene Summary: This gene belongs to subfamily P of the forkhead box (FOX) transcription factor family. Forkhead box transcription factors play important roles in the regulation of tissue- and cell type-specific gene transcription during both development and adulthood. Many members of the forkhead box gene family, including members of subfamily P, have roles in mammalian oncogenesis. This gene may play a role in the development of tumors of the kidney and larynx. Alternative splicing of this gene produces multiple transcript variants, some encoding different isoforms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG215338