

Product datasheet for **RG215462**

OR6C65 (NM_001005518) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	OR6C65
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PSI00010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG215462 representing NM_001005518
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGCCAAATATGACATCAATTAGAGAATTCATTCTTCTGGGATTTACAGATAACCCAGAGTTACAAGTTG
TGATATTCCTTTATGTTGATCACATACTATTGAGTGTAAGTGAAACATGATCATCATTATGTTAAC
ATTGTCAAATATTCATTTGAAAACCTCTATGTATTTCTTCCTCAGGAATTTCTTTTCTTAGAAATTTCA
TTTACGACAGTCTTCATTCTAGATTTCTGATCAACATTGCTACAGGAGACACAACCATTTCTATAATG
CTTCCATGGCCCAAGTATTTTTTAATCTTTTGGGATCAACAGAATTTTTCTTCTGGCTGTCATGTC
TTATGATCGCTATGTGGCTATCTGCAAACCTCTACATTATAACAACCATCATGAGTAACAAAGTCTGTAAT
TGGCTTGTAATCAGCTCCTGGCTGGCTGGTTTTCTCATTATTTTTCCCCCGTGATTATGGGCTCCAAC
TGGATTTTTGTGACTCCAGCACTATTGACCATTTCATCTGTGATTCTCCCTATGCTGCTGATTGCTTG
CACAGACACAGTTTCTAGAGCTTATGGCATTCTTGCTAGCAGTATTCACACTCATGGTAACCTTTGGCC
TTGGTGGTTCTCTCTACACTCATCCTTAAACAATTCTGAAGATTCCTCTGCCAGCAAAGGAAAA
AGGCTTTTTCAACTGTTCCTCCCATATGATTGTGTTTCTGTTTCTTATGGGAGTTGCATTTTATGTG
TGTA AAAACATCTGAAAAGAAGGTATGGCTTTGAGCAAAGGTGTAGCAGTGCTTAATACCTCTGTTGCT
CCTATGTTGAATCCCTTTATTTATACCTTAAGAAACCAGCAGGTGAAACAGGCCCTTAGGGAATTCACCA
AAAAATATTATCATTGAACAAGCAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence:

>RG215462 representing NM_001005518

Red=Cloning site Green=Tags(s)

MPNMTSIREFILLGFTDNPELQVVIFFMLITYLLSVSGNMIIIMLTLSNIHLKTPMYFFLRNFSFLEIS
FTTVFIPRFLINIATGDTTISYNASMAQVFFLLLGSTEFFLLAVMSYDRYVAICKPLHYTTIMSNKVCN
WLVISSWLAGFLIIFPPVIMGLQDFCDSSTIDHFCIDSSPMLLIACDTQFLELMAFLLAVFTLMVTLA
LVVLSYTLILKILKIPSAQQRKKAFTSCSSHMIVSVSYGSCIFMCVKTSAKEGMALSKGVAVLNLSVA
PMLNPFITYTLRNQYKQALREFTKILSLNKQ


TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



CTATAGGGCGGCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGATCGCCGGCGCCAGATCT

EcoRI BamHI KpnI RBS Kozac Consensus SgfI AscI

CAAGCTTAAGCTAGCTAGCGGACCG ACG CGT ACG CGG CCG CTC GAG ATG GAG AGC GAC ---

HindIII NheI RsrII MluI NotI XhoI GFP Tag

T R T R P L E M E S D - - -

--- GAA GAA AGA GTT TAA ACGGCCGGCCGCGGAGCT

- - E E R V Stop

ACCN:

NM_001005518

ORF Size:

936 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:


This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001005518.1, NP_001005518.1</u>
RefSeq Size:	939 bp
RefSeq ORF:	939 bp
Locus ID:	403282
UniProt ID:	<u>A6NJZ3</u>
Cytogenetics:	12q13.2
Protein Families:	Transmembrane
Protein Pathways:	Olfactory transduction
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG215462