

Product datasheet for **RG215478**

CEACAM18 (NM_001080405) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	CEACAM18
Synonyms:	2010110O04Rik; carcinoembryonic antigen-related cell adhesion molecule 1; carcinoembryonic antigen-related cell adhesion molecule 18; CEA-related cell adhesion molecule 1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG215478 representing NM_001080405
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGGAGGTGTCTGGCTGCCCCACCCCTGCTGGTCAGAGTGCACCCTGGGGAGGAAAGGCCGCCCTGAAGG
CTGGGAGACGAGACCGGCAGCCTCTGTGCCAGGAGACACAGGTCTTGGAGGGAGCAGAGGAGGTGGCTGT
GTCTCTGGAGGGACCCCTCCTCTGGAGGACCCAGGCAGCTCATGGACCTTTCAGACCCAGATGGAGC
CTGTGGAGGAGGGTCTTCTCATGGCCAGTCTGCTGGCCTGTGGGATCTGCCAGGCCTCTGGCCAAATCT
TCATCACCCAAACCCTGGGGATCAAGGGATATCGGACTGTCTGGCCCTGGATAAGGTCCCTGAGGATGT
TCAGGAATACAGCTGGTACTGGGGTGCAAACGACAGCGCAGGAAACATGATTATCAGCCACAACCCGCC
AGTGCCAGCAGCCTGGGCCATGTACACTGGCAGGAGAGAGTGAACAGAGAAGGCAGCCTGTTGATCA
GGCCGACTGCATTAAATGACACGGGAACTACACTGTTCCGGTGGTTCAGGCAATGAGACCCAAAGAGC
AACCGGCTGGCTGGAGTTCTAGAGTTGGGAAGCAATCTGGCATCTCCGTCAATGCCAGCTCCCTGGTG
GAGAACATGGATTCTGTGGCTGCTGACTGCCTCACAATGTACCAACATCACGTGGTATGTGAATGATG
TGCCTACCTCTAGTAGTGACCGGATGACAATTTCCCCAGACGGCAAGACCTCGTCATCCTCAGGGTCAG
CCGCTATGACAGACAATTTCAGTGCATGATAGAGAGTTCCAGAGATCTTTCAGAGAAGTGAACGCATC
TCTCTGACTGTGGCCTATGGGCCGACTATGTGCTGCTGAGGAGCAATCCTGATGATTTCACCGCATTG
TGACAGCTGAGATCGGCTCCCAAGTGGAATGGAGTGTATCTGCTATTCCCTTCCTGGATCTCAAGTACCA
CTGGATCCACAATGGCTCCCTCCTGAATTCTCAGATGCAAAGATGAACCTCTCGAGTCTTGCCTGGGAG
CAGATGGGCCGTTACCGATGCACTGTGGAGAACCCGTGACACAGCTGATCATGTACATGGACGTCAGGA
TCCAGGCCCCCATGAGTGCAGCAGCTCCCTCCAGGCTCATGCTTGCACATCTCCTGCCTCCATGCC
CTGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG215478 representing NM_001080405
Red=Cloning site Green=Tags(s)

MEVSGCPTPAGQSAPWGGAALKAGRDRQLCQETQVLEGAEVAVSLEGPLLLEDPRQLMDLSRPRWS
LWRRVFLMASLLACGICQASGQIFITQTLGIKGYRTVVALDKVPEDVQEYSWYWGANDSAGNMIIISHKPP
SAQQPGPMYTGREVRNREGSLLRPTALNDTGNYYTRVVAGNETQRATGWLEVLGLSNLGISVNASSLV
ENMDSVAADCLTNVTNITWYVNDVPTSSSDRMTISPDGKTLVILRVSRYDRTIQCMIESFPEIFORSERI
SLTVAYGPDYVLLRSNPDFNGIVTAIEGSCQVEMECISFYSLDLKYHIHWIHGSLNLFSDAKMNLSSLAWE
OMGRVRYCTVENPTOLIMYDNRTOAPHEGSSSPPGSCFAHLPASMPCC

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001080405

ORF Size: 1194 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001080405.1](#), [NP_001073874.1](#)


RefSeq Size: 1197 bp

RefSeq ORF: 1196 bp

Locus ID: 729767

Cytogenetics: 19q13.41

Product images:



Circular map for RG215478