

Product datasheet for **RG215482**

PSMF1 (NM_178578) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PSMF1 (NM_178578) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PSMF1
Synonyms:	PI31
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG215482 representing NM_178578 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGGCCTGGAGGTACTGTTTCGCATCGGCAGCGCCGCCATCACCTGCAGGCAGGACGCGCTCGTCT
GCTTCTTGCATTGGGAAGTGGTGACACACGGTTACTGCGGCTTGGGTGTCGGTGACCAGCCGGTCCCAA
TGATAAGAAGTCAGAAGTCTGCCAGTGGGTGGAACAACAATAAAGACCTGTATGTCTCCGGTATGAG
TATAAGGATGGTCCAGAAAGCTCCTTGTGAAAGCCATCACCGTGGAGAGCAGCATGATCCTCAATGTGC
TGAATATGGCTCACAGCAAGTGGCAGACTTGACCCTGAACTTGGATGATTATATCGATGCAGAACCT
GGGTGACTTCCACAGGACCTACAAGAACAGTGAGGAGCTTCGGTCTCGTATTGTGTCTGGAATCATCACA
CCTATCCATGAGCAGTGGGAAAAGGCTAATGTAAGCAGTCCCCACCGGGAGTTCCCCCTGCTACCGCCA
GAGAGGTGGACCACTCCGGATTCTCCACACCACCCACACACCAGTCGGCAGCCTCCCTGGTGTGATCC
CCTGGGCCGTTTGTTCGGGGGAGAAGACTTAGACCCTTTTGGGCCTCGGAGAGGTGGCATGATTGTG
GATCCCCTGAGATCTGGCTTCCAAGAGCACTTATTGACCCTTCTCAGGCCTCCCGAACCGACTTCCTC
CAGGCGCTGTGCCCCAGGAGCTCGCTTGGACCCTTTGGACCATTGGGACCAGCCACCCGGACCTAA
CCAGACCATCTCCCCCGCCGGCTACGATGACATGTACCTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG215482 representing NM_178578
 Red=Cloning site Green=Tags(s)

MAGLEVLFAAAPAITCRQDALVCFLHWEVVTHGYCGLGVDQPGPNDKSELLPAGWNNKDLVYLRYE
 YKDGSRKLLVKAITVESSMILNVLEYGSQQVADLTNLDDYIDAEHLGDFHRTYKNSEELRSRIVSGIIT
 PIHEQWEKANVSSPHREFPPATAREVDPLRIPPHHPHTSRQPPWCDPLGPFVVGGEDLDPFGPRRGGMIV
 DPLRSGFPRALIDPSSGLPNRLPPGAVPPGARFDPFPGPIGTSPPGPNPDLPPPGYDDMYL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_178578

ORF Size: 813 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_178578.1](#), [NP_848693.1](#)

RefSeq Size: 3204 bp

RefSeq ORF: 816 bp

Locus ID: 9491


UniProt ID: [Q92530](#)

Cytogenetics: 20p13

Protein Pathways: Proteasome

Gene Summary: The 26S proteasome is a multicatalytic proteinase complex with a highly ordered structure composed of 2 complexes, a 20S core and a 19S regulator. The 20S core is composed of 4 rings of 28 non-identical subunits; 2 rings are composed of 7 alpha subunits and 2 rings are composed of 7 beta subunits. The 19S regulator is composed of a base, which contains 6 ATPase subunits and 2 non-ATPase subunits, and a lid, which contains up to 10 non-ATPase subunits. Proteasomes are distributed throughout eukaryotic cells at a high concentration and cleave peptides in an ATP/ubiquitin-dependent process in a non-lysosomal pathway. An essential function of a modified proteasome, the immunoproteasome, is the processing of class I MHC peptides. This gene encodes a protein that inhibits the activation of the proteasome by the 11S and 19S regulators. Alternative transcript variants have been identified for this gene. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG215482