

Product datasheet for **RG215511**

CCDC81 (NM_021827) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC81 (NM_021827) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CCDC81
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG215511 representing NM_021827
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**C

ATGTTGGATACGATCGCCCGTGCCCTGCAGGACCTGGGCAGGCAGGTGCTGCCACTCTGCCCTCGCTGA
 GCCAGGAGGAAGTCTCTATTATCTGGGGGAATGTATCAGAATTTGTGAGACGGCAGTTAACCTGCACAA
 GGGGGTTCAGATTCCAGCATTTGGAACCTTCACCTTCATAAGACAAAAGCTTGAGGTGGGAAACAAACAA
 TTTATCTTAATCCAGAGGCTGTGTTTATCATGGTGGGAAGCTAGTGCAGATCCATGGACTCAAACAA
 ACAAAGGCTGGCACTGTGGACTCGGTGTTGTCTAGCAGAGAGGCTTGAGGAAGTGGCCAGCAGTGT
 GCTTGCCTTTCCAAGGATTGAGCTCAAGGAGATGGAAAACAACTGCCTATGGAGACCCTGTAGAAGAA
 TGTGGAGAGAATAGAGAAAGGAAGTGCAAGTTAAAGACCAGTCAGACAAAGAAGAAGGCACCAGAGATA
 TCTCATACCCAAAAGACTTCGAGATAGACAAGCTTTGTTCCCTGCCAAAGTGACAAATGTCAGCTTGCT
 GGGAAAGTTTGAACGAAGTGAGAGTGGTGGGAAGATTATGACCCCTGAAAGCTTATCATATCCAAGTTGT
 CTGAAACACGACAGTGAGATGAAGCCCCAACATCTCCAGCTTGCCAGGATCATAAACAGGCAGGACAGG
 AAATGTGCTATGTATGTTTGCAACGAGCACAACGAAATTCCTGTGTACTACAGTGAGGAAGGAGGAG
 AGAGATAGAAGATGAGAGACTCATACAGCAGTATCAGATGTTAAAGGATCAGGAGGCTCTCTTCAGACAC
 CAGATGAAAAGTCTGGCTACTAGAGAACAGAATCAGAAAAATGCTGCCTATAATCTTGGAGTTGCTGAAG
 CTATAAGAAACCAAGAATGAGAAACCGGAATTTTATAAATCCTTCCTATTTGACAAACGGCCACTCAG
 TCCTGCGCTTAATGCTCTTAAGCAAGAGGAATATTCGGAGTCTCCTGAAACAAATGGATAACAGACAG
 GAAAACGAAATAAAGCAAAGACAATACAGAGAGTTGATGGACCGCTGGAACAAGTGCAACTCACAGAGG
 AACTTGCTGCGCAAAGAGCGAAATTTTAAAAGATAAGATGGAAGAAACACAGTGTTACAAGAGAGCTTT
 GGATGCACAGATAAAGAACAACCTCTCGGCTGCCCCCTTTGAGCCAGACTCCTCTGAGCCCATCTTT
 GGTAAGAATGAGGTGAACTGATGGTGGAAAAGCAAAAGCGAGAACAAATACATGAAACACCACTGG
 AGGCAGCTGCTAACCACAAGAGGAAAGCCATCCTGCATCAACTAGTGGACCAGAGCGGGATTTGCAAA
 GCTTCAGAGGACACAAAGAGAGCACTTGGCAGACAGAACCCTGAGCTGGAGCGAGTAAATAGAGTCAAC
 CAATGCTTACAGGAGGACTGGGAAAGGAGTGTGCGATGAAGAAGCAGCGAGACCTGGAGGACAAGGCTT
 TTGAACGGGCTTCAGACAAGCTGTTTCTCTAGACCAGTGTGAGAAGTATCGGCGCTGCAAGCAATGCCA
 GAGGCGCACCTCCAACGTGGGCGAGAGCAACCTGTGGCCCTGAACAAGTTCCTGCCTGGCTCCCGGTTG
 CTTGTG

ACGCGTACGCGGCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG215511 representing NM_021827
 Red=Cloning site Green=Tags(s)


MLDTIARALQDLGRQVLPTLPSLSQEEVSIWGNVSEFVRRQLTLHKGVIQIPAFGFTTFIRQKLEVGNNK
 FILIQRPVFIWVGLVQIHGLKQNKRPVTVDSVLSREALRKWPSSVLAFPRIELKEMENKLPMETLVEE
 CGENRERKCKLKDQSDKEEGTRDISPKRLRDRQALFPAKVTVNSLLEKFERSESQKIMTPESLSYPSC
 LKHDSEMPQTPACQDHNKAGQEMCYVCLQRAQRNSLLYYSEERRREIEDERLIQQYQMLKDQEALEFRH
 QMKSLATREQNQNAAAYNLGVAEAIKRNHNEKPEFYKSLFDKRPLSPALNALKQEEYSRLLKQMDNRQ
 ENEIKQRQYRELMDRLEQVQLTEELAAQRAKFLKDKMEETQCYKRALDAQIKNKPSRLPPFEPDSSEPIF
 GKNEGELMVEKQKREQNYMKHQLEAAANHKRAILHQLVDQRRDLQMLQRTQREHLADRTAELERVNRVN
 QCLQEDWERSAAMKKQRDLEDKAFERASDKLFLLDQCEKYRRCKQCQRRTSNVGESNLWPLNKFPGSRL
 LV

TRTRPLE – GFP Tag – V

Restriction Sites:


Sgfl-MluI

Cloning Scheme:



ACCN:	NM_021827
ORF Size:	1686 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_021827.3
RefSeq Size:	2513 bp
RefSeq ORF:	1689 bp
Locus ID:	60494
UniProt ID:	Q6ZN84
Cytogenetics:	11q14.2

Product images:



Circular map for RG215511