

Product datasheet for **RG215564**

TSPYL4 (NM_021648) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	TSPYL4
Synonyms:	dJ48613.2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG215564 representing NM_021648
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGAGCGGCCTGGATGGGGCAACAAGCTCCCTCTCGCCAAACCGCGGCCTGGCTGCTCCCGACCATG
CCTCAGGAGATCCGGACCGAGACCAAGGCTCCGTGAAGAAACCGAGGCGACACAGGTGATGGC
GAACACAGGTGGGGGACGCTGGAGACCGTTGCGGAGGGGGTGCATCCAGGATCCTGTCGACTGTGGC
CCCGCGCTCCGCGTCCAGTTGCCGGGAGTCGCGCGGTGCAGCGACCAAGCCGGGCAGGAGGATGCTC
CACCTTCTACGAAAGGTCTGGAAGCAGCCTCTGCCGCCGAGGCTGCTGACAGCAGCCAGAAAAATGGCTG
TCAGCTTGAGAGCCCCGTGGCCCTGCTGGGCAGAAGGCTCTAGAAGCCTGTGGCGCAGGGGGCTTGGGG
TCTCAGATGATACCGGGAAGAAGGCCAAGGAAGTGACGACTAAAAACGCGCCATCTCGGCAGCAGTGG
AAAAGGAGGAGAAGCAGGGGCGCGATGGAGGAAAAGAAGTAGTGAGAAGGAAAAAAGGTGGCAGG
AGGGGTGAAAGAGGAGACACGGCCAGGGCCCCGAAGATCAATAACTGCATGGACTCACTGGAGGCCATC
GATCAAGAGTTGTCAAACGTAAATGCCAGGCTGACAGGGCCTTCCTTCAGCTTGAGCGCAAGTTTGGCC
GCATGCGAAGGCTCCACATGCAGCGCAGAAGTTTCAATATCCAGAATATCCAGGTTTCTGGGTTACTGC
CTTTCGAAACACCCAGCTGTACCTATGATCAGTGCCCAAGATGAAGACATGCTGAGGTACATGATC
AATTTGGAGGTGGAGGAGCTTAAACACCCAGAGCAGGCTGCAATTCAAGTTCACTTTTCAGGGCAACC
CCTACTTCCGAAATGAGGGGCTTGTCAAGGAATATGAACGCAGATCCTCTGGCCGGGTGGTGTCTTTTC
CACTCCAATCCGCTGGCACCAGGCAAGACCCAGGCTCATATCCACAGAAACGGGAAGGGAACACT
ATCCCTAGTTTCTCAACTGGTTTTAGACCACAGCCTTCTAGAATTCGACAGAATTGCAGAGATTATCA
AAGGAGAACTGTGGCCCAATCCCTACAATACTACCTGATGGGTGAAGGGCCCCGTAGAGGAATTCGAGG
CCCACCAAGGCAGCCAGTGGAGAGCGCCAGATCCTTCAGGTTCCAGTCTGGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG215564 representing NM_021648
Red=Cloning site Green=Tags(s)

MSGLDGGNKLPLAQTTGGLAAPDHASGDPDRDQCQGLREETEATQVMANTGGGSLETVAEGGASQDPVDCG
PALRVPVAGSRGGAATKAGQEDAPPSTKGLEAASAAEADSSQKNGCQLGEPRGPAGQKALEACGAGGLG
SQMIPGKKAKEVTTKKRAISAAVEKEGEAGAAMEEKKVVQKEKKVAGGVKEETRPRAPKINNCDMSLEAI
DQELSNVNAQADRAFLQLERKFGMRRLHMQRSSFIIQNIPGFVWTAFRNHPQLSPMISGQDEDMRLRYMI
NLEVEELKHPRAGCKFKFIFQGNPYFRNEGLVKEYERRSSGRVVSLSLTPIRWHRGQDPQAHIRNREGNT
IPSFFNWFSHSLLEFDRIAIEIIKGLWPNPLQYYLMGEGPRRGIRGPPRPQPVESARSFRFQSG

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_021648

ORF Size: 1242 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_021648.5

RefSeq Size: 4144 bp


RefSeq ORF: 1245 bp

Locus ID: 23270

UniProt ID: Q9UJ04

Cytogenetics: 6q22.1

Product images:



Circular map for RG215564