

## Product datasheet for **RG215704**

### PPM1B (NM\_001033557) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PPM1B (NM_001033557) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PPM1B
Synonyms:	PP2C-beta; PP2C-beta-X; PP2CB; PP2CBETA; PPC2BETAX
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG215704 representing NM_001033557 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGGTGCATTTTTGGATAAACCCAAAACCTGAAAAACATAATGCTCATGGTGTGGGAATGGTTTACGTT  
ATGGCCTGAGCAGCATGCAAGGATGGAGAGTGGAAATGGAAGATGCACACACAGCTGTTGTAGGTATTCC  
TCACGGCTTGAAGACTGGTCATTTTTGTCAGTTTATGATGGTCATGCTGGATCCCGAGTGGCAAATTAC  
TGCTCAACACATTTATTAGAACACATCACTACTAACGAAGACTTTAGGGCAGCTGGAAAATCAGGATCTG  
CTCTTGAGCTTTCAGTGGAAAATGTTAAGAATGGTATCAGAACTGGATTTTTGAAAATTGATGAATACAT  
GGTAACCTTTTCAGACCTCAGAAACGGGATGGACAGGAGTGGTTCAACTGCAAGTGGGAGTTATGATTTCA  
CCTAAGCATATCTACTTTATCAACTGTGGTGATTCACGTGCTGTTCTGTATAGGAATGGACAAGTCTGCT  
TTTCTACCCAGGATCACAAACCTTGCAATCCAAGGGAAAAGGAGCGAATCCAAAATGCAGGAGGCAGCGT  
GATGATACAACGTGTTAATGGTTCATTAGCAGTATCTCGTCTGGGGACTATGATTACAAGTGTGTT  
GATGGCAAGGGCCCAACAGAACAACCTGTTTCTCCAGAGCCTGAGGTTTATGAAATTTAAGAGCAGAAG  
AGGATGAATTTATCATCTGGCTTGTGATGGGATCTGGGATGTTATGAGTAATGAGGAGCTCTGTGAATA  
TGTTAAATCTAGGCTTGAGGTATCTGATGACCTGGAAAATGTGTGCAATTGGGTAGTGGACACTTGTTA  
CACAAGGGAAGTCGAGATAACATGAGTATTGTACTAGTTTGCTTTTCAAATGCTCCCAAGGTCTCAGATG  
AAGCGGTGAAAAAGATTACAGATTGGATAAGCACTTGGAAATCACGGGTTGAAGAGATTATGGAGAAGTC  
TGGCGAGGAAGGAATGCCTGATCTTGCCCATGTATGCGCATCTTGCTGTCAGAAAAATATCCCAAATTTG  
CCTCCTGGGGAGGTCTTGCTGGCAAGCGTAATGTTATTGAAGCTGTTTATAGTAGACTGAATCCACATA  
GAGAAAGTATGGGCAGAAA

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online >](#)

Protein Sequence: >RG215704 representing NM\_001033557  
 Red=Cloning site Green=Tags(s)

MGAFLDKPKTEKHNAHGAGNGLRYGLSSMQGWRVEMEDAHTAVVGIPHGLEDWSFFAVYDGHAGSRVANY  
 CSTHLLLEHITTNEDFRAAGKSGSALELSVENVKNGIRTGFLKIDEYMRNFSDLRNGMDRSGSTAVGMIS  
 PKHIYFINCGDSRAVLYRNGQVCFSTQDHKPCNPREKERIQNAGGSVMIQRVNGSLAVSRALGDYDYKCV  
 DGKGPTEQLVSPPEVYEILRAEEDEFIILACDGIWDVMSNEELCEYVKSRLVSDLENVNWNVVDTCLE  
 HKGSRDNMSIVLVCFSNAPKVSDEAVKDSLELDKHLSESRVEEIMEKSGEEMGMPDLAHVMRILSAENIPNL  
 PPGGLLAGKRNVIEAVYSRLNPHRESDGQK


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001033557

<b>ORF Size:</b>	1140 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001033557.3</a>
<b>RefSeq Size:</b>	3056 bp
<b>RefSeq ORF:</b>	1143 bp
<b>Locus ID:</b>	5495
<b>UniProt ID:</b>	<a href="#">O75688</a>
<b>Cytogenetics:</b>	2p21
<b>Protein Families:</b>	Druggable Genome, Phosphatase, Stem cell - Pluripotency
<b>Protein Pathways:</b>	MAPK signaling pathway
<b>Gene Summary:</b>	The protein encoded by this gene is a member of the PP2C family of Ser/Thr protein phosphatases. PP2C family members are known to be negative regulators of cell stress response pathways. This phosphatase has been shown to dephosphorylate cyclin-dependent kinases (CDKs), and thus may be involved in cell cycle control. Overexpression of this phosphatase is reported to cause cell-growth arrest or cell death. Alternative splicing results in multiple transcript variants encoding different isoforms. Additional transcript variants have been described, but currently do not represent full-length sequences. [provided by RefSeq, Jul 2008]