

Product datasheet for **RG215711**

LIN54 (NM_194282) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	LIN54
Synonyms:	CXCDC1; JC8.6; MIP120; TCX1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



ORF Nucleotide Sequence: >RG215711 representing NM_194282
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGAGGTGGTGCCAGCTGAGGTGAGTAGTTTGCTTCCAGAGGAAATAATGGACACTGGTATAACTTTAG
 TGGATGATGATAGTATTGAGGCTGTTATTGTTTCATCCCAATTCCCATGGAGACAGAACTGGAAGAAAT
 TGTCACATAAAATCTACTGGTGACTCTACAGCCACGCCCATTTCACGGAACCAATCACAGTGTACAGT
 AACCACACTAACCAAGTTCAGTGAATACCACAATTACTAAGCAGATTCTAATACCACAGTGAAACCAG
 CTTTTCCAAGTGGCCTTCAAAAACCTTGGTGCTCAGACTCCTGTGACTATATCAGCCAATCAGATTATTTT
 AAACAAAGTATCAGACATCTGATCTTAAACTTGGCAATCAGACCCTTAAACCAGATGGACAGAAGTTA
 ATTTTAAACACTTTGGGCAAGTCTGGTTCACCAATTGTTTTAGCACTACCCATAGCCAACCTACCCAGG
 CTCAGAAAGTTACAACCTCAGGCCAGTCAGGAGATGCTAAGTTACCACCGCAGCAAATTAAGTAGTTAC
 CATTGGAGGGAGCCAGAGGTGAAACCTGTCATTGGTGTCTCAGCATTGACCCAGGAAGTCAACTGATT
 AATACTACAACCTCAGCCCTCTGTGTACAGACCCAACAGTTAAAAACAGTACAGATTGCTAAGAAGCCTC
 GAACGCCAACCTCTGGTCCAGTAATCACGAAGCTGATCTTTGCAAAACCAATTAATAGTAAAGCAGTTAC
 AGGACAGACAACCTCAAGTTTACCACCAAGTCATTGCAGGTAGGGTCTTTTACAGTCTACTCCCGGAAC
 CCATCAAAGACCATAACAATATCTGAAAGTGGTGTATTGGATCACTTTAAATTTTACAACACAGACAC
 CAAATAAAATAGCCATCTCACCTTTGAAATCGCCAAATAAGGCAGTGAAATCAACTGTGCAGACCATCAC
 GTTGGAGGAGTGAGCAGATCAGGTTTAAAGACAATTATTCCTCTGGCAACTGCTCCCAATGTCCAGCAG
 ATTCAGTGCCTGGAAGCAAGTTTCATTATGTCCGACTTGTACTGCCACATCAGCCAGTAGCTCAACCC
 AGCCAGTTAGTCAGAAATCCAGTACAAACACTCAGCCTCTTCAGCAAGCAAAGCCAGTGGTTGTTAATAC
 AACCCAGTGCGGATGTGAGTTCCAATTGTCTCAGCTCAGGCTGTCAAACAAGTTGTTCCAAAACCAATC
 AATCCAACCTCACAAATAGTAACTACTAGCCAGCCACAGCAACGGCTTATCATGCCCCGCACACCACTGC
 CACAGATCCAGCCCAACCTCACTAACCTGCCACCAGGCACTGTCTGGCACCAGCTCCGGGAACAGGGAA
 TGTGGGTTATGCACTGCTTCCAGCTCAGTATGTTACTCAGCTACAGCAGTCTTCATATGTATCAATAGCA
 AGCAACTCTACCTTTACTGGAACATCTGGTATCCAGACCCAGGCACGGCTTCCATTCAATGGCATAATCC
 CATCAGAGTCGGCCAGTCGGCCCCGAAAGCCCTGTAATTGTACAAAATCACTGTGTTTGAAATTGTATTG
 TGATTGCTTTGCAAAATGGTGAATTTTGCAACAACCTGCAATTGTACTAATTGTTACAACAATTTGGAACAT
 GAAATGAAAGGCAAAAAGCAATAAAGGCATGCCTTGACAGAAATCCAGAAGCCTTTAAGCCTAAGATAG
 GGAAAGGAAAGGAGGGAGAATCTGATCGACGTCATAGCAAAGGGTGAATTGCAAACGATCAGGATGTCT
 TAAAACTACTGTGAATGCTATGAGGCAAAAATAATGTGTTCTCAATATGCAAAATGTATTGGCTGTAAG
 AATTTTGAAGAAAGCCCGAAAGGAAGACATTGATGCATTTGGCAGATGCAGCTGAAGTAAGGGTACAGC
 AACAAACAGCAGCCAAGACGAAGTTATCCTCTCAAATTTTCAAGTTGCTTACTAGGCCAACACCACTTT
 AAATAGTGGAGGCGGAAAATTGCCATTTACATTTGTAACCTAAGGAAGTAGCTGAAGCCACATGTAATTGC
 CTCCTTGCCAGGCAGAGCAGGCAGACAAGAAGGGAATCAAGGCAGCAGCGGAACGGATGATACTTG
 AGGAATTCGACGATGTTTGATGAGTGTCACTCTGCAGGAAAGGCAAAAAGTGCCTTGTGCCAT
 GAATTGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG215711 representing NM_194282

Red=Cloning site Green=Tags(s)


MEVVP AEVSSLLPEEIMDTGITLVDDDSIEAVIVSSPIPMETELEEIVNINSTGDSTATPISTEPITVYS
NHTNQVAVNTTITKADSNNTTKPAFPSPGLQKLGAQTPVTISANQIILNKVSQTSCLKLGNQTLKPDGQKL
ILTTLGKSGSPIVLALPHSQLPQAQKVTTQAQSGDAKLPPQQIKVVTIGRPEVKPVIGVSALTPGSQLI
NTTTQPSVLQTTQQLKTVQIAKKPRTPTSGPVITKLIFAKPINSKAVTGQTTQVSPPIAGRVLSQSTPGT
PSKTIITISESGVIGSTLNSTTQTPNKIAISPLKSPNKAVKSTVQTIITVGGVSTSQFKTIIPLATAPNVQQ
IQVPGSKFHYVRLVTATSASSSTQPVSNPSTNTQPLQQAQKPVVNTTPVRMSVPIVSAQAVKQVVPKPI
NPTSQIVTTSQPPQRLIMPATPLPQIQPNLTNLPPGTVLAPAGTGNVGYAVLPAQYVTQLQSSYVSIA
SNSTFTGTSGIQTQARLPFNGIIPSESASRPKPCNCTKSLCLKLYCDFANGFECNNCNCNTNCYNNLEH
ENERQKAIKACLDNRPEAFKPKIGKKEGESDRRHSGKCNCKRSGCLKNYCECYEAKIMCSSICKIGCK
NFEESPERKTLMLADAAEVRVQQQTAATKLSQISDLLTRPTPALNSGGGLPFTFVTKEVAEATCNC
LLAQAEQADKKGSKAAAERMILEEFGRCLMSVINSAGKASAPCAMNC

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_194282

ORF Size:

2247 bp

OTI Disclaimer:


The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_194282.1 , NP_919258.1
RefSeq Size:	5772 bp
RefSeq ORF:	2250 bp
Locus ID:	132660
UniProt ID:	Q6MZIP
Cytogenetics:	4q21.22
Gene Summary:	LIN54 is a component of the LIN, or DREAM, complex, an essential regulator of cell cycle genes (Schmit et al., 2009 [PubMed 19725879]).[supplied by OMIM, Dec 2010]

Product images:



Circular map for RG215711