

Product datasheet for **RG215741**

PRELID2 (NM_182960) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: PRELID2 (NM_182960) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: PRELID2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG215741 representing NM_182960
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC**CGCATCGCC**

ATGGGGGTCTCGGTGGATGTGCACCAAGGTGTACAAGTACCCCTTCGAGCAGGTGGTCGCCAGCTTTCTCC
 GAAAGTACCCCAACCCCATGGATAAAAAATGTCATCTCAGTAAAAATCATGGAGGAAAAAGAGATGAATC
 AACAGGGGTCACTACAGAAAGAGGATTGCAATCTGTGAGAACGTGGTTCAGAAATTTAAGGAAGTCC
 TTAAGCACTCTGGTGATTCTTTGTTGGAAAAAGGTGAGCATTTTGAAAGTACCTAATATCCAATTAGAAG
 AGGAGTCATGGCTCAATCCTCGGAAAGAAACATGGCCATACGGAGTCACTGCCTTACGTGGACACAGTA
 TGCATCCATGAAGGAAGAGTCTGTCTTCCGGGAAAGTATGGAAAACCCAAATTGGACAGAGTTCATTCAA
 AGAGGCAGGATTTCAATCACAGGGGTTGGATTTCTCAACTGTGTTTTAGAACTTTGCCAGCACATTCT
 TACGACAGGGAGCCAGAAAGGAATTAGAATCATGGAGATGCTGCTAAAGGAACAGTGTGGTGCCCCCTT
 AGCTGAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG215741 representing NM_182960
 Red=Cloning site Green=Tags(s)

MGVSVDVHQVYKYPFEQVVASFLRKYPNPMKDNVISVKIMEEKRDSTGVIYRKRIAICQNVVPEILRKS
 LSTLVILCWKKVSILKVPNIQLEESWLNPRERNMAIRSHCLTWTQYASMKEESVFRESMENPNWTEFIQ
 RGRISITGVGFLNCVLETFASTFLRQGAQKIRIMEMLLKEQCGAPLAE

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: NM_182960

ORF Size: 567 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_182960.4](#)

RefSeq Size: 2140 bp


RefSeq ORF: 570 bp

Locus ID: 153768

UniProt ID: [Q8N945](#)

Cytogenetics: 5q32

Product images:



Circular map for RG215741