

Product datasheet for **RG215748**

PTPN20B (PTPN20) (NM_001042357) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PTPN20B (PTPN20) (NM_001042357) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PTPN20
Synonyms:	bA42B19.1; bA142I17.1; CT126; PTPN20A; PTPN20B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG215748 representing NM_001042357 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTTCACCTAGGGACTTTAGAGCAGAGCCTGTAACGATTATGAGGGAAATGACTCTGAAGCAGAAG
ACTTGAATTTAGGGAGACTTTGCCTTCATCAAGTCAGGAAAACACACCTAGATCAAAGGTTTTTAAAA
TAAAGTTAATTCAGAGAAGGTAACCTTCTCTCGGAATTTCCACATAATGATTATGAGGATGTTTT
GAAGAGCCTTCAGAAAGTGGCAGTGATCCAGCATGTGGACAGCCAGAGGCCCTTCAGAAGAGACAGGT
GGAGCAGTGAGGATGAGGAGGCTGCAGGGCCATCACAGGCTCTCTCCCTCTACTTTCTGATACGCGCAA
AATTGTTTCTGAAGGAGAACTAGATCAGTTGGCTCAGATTCGGCCATTAATATTCATTTTCATGAGCAG
ACAGCCATCAAGGATTGTTTAAAAATCCTTGAGAAAAAACAGCAGCGTATGATATCATGCAAGAAATTTA
TGGCTTTAGAAGTAAAGAACTGCTGCTGGTGAAGTCAACTCTGGGAATCAACCAAGCAACAGAGAAAAAA
TAGATACCGAGATATTCTTCCATATGATTCAACACGCGTTCCTCTTGAAAAAGCAAGGACTACATCAAT
GCTAGTTATATTAGAATAGTCAATTGTGGAGAAGAGTATTTTTATATCGCTACTCAAGGACCACTGCTGA
GCACCATAGATGACTTTTGGCAAATGGTGTGGAAAAATAATCAAATGTTATTGCCATGATAACCAAGAGA
GATAGAAGGTGGAATTATCAAATGCTACCATTACTGGCCATTTCTCTGAAGAAGCCATTGGAATTGAAA
CACTTCCGTGATTCTCTGGAGAAGTACCAGATACTTCAATATTTTCATTCGAAATGTTTCAAGTTGTGG
AGAAGTCCACGGGAAGTCACTCTGTAACACAGTTGCAAGTTCACCAAGTGGCCAGACCATGGCACTCC
TGCTCAGCAGATAGCTTCATAAAATATATTCGTTATGCAAGGAAGGCCACCTTACAGGACCCATGGTT
GTTCACTGCAGTGCCGGCATAGGCCGGACAGGGGTGTTCTATGTGTGGATGTCGTGTTCTGTGCCATCG
TAAAGAAGTGTTCATTCAACATCATGGATATAGTGGCCAAATGAGAGAACAACGTTCTGGCATGGTTCA
AACGAAGGAGCAGTATCACTTTTGTACGATATTTGCTTGAAGTTCTTCGAAACTTCTGACTTTGGAT

AGCGGACCGACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG215748 representing NM_001042357
 Red=Cloning site Green=Tags(s)

MSSPRDFRAEPVNDYEGNDSEAEDLNFRETLPSSSQENTPRSKVFENKVNSEKVKLSLRNFPNDYEDVF
 EEPSESGSDPSMWTARGPFRDRWSEDEEAAGPSQALSPLLSDTRKIVSEGELDQLAQIRPLIFNFHEQ
 TAIKDCLKILEEKTAAYDIMQEFMALELKNLPGEFNSGNQPSNREKNRYRDILPYDSTRVPLGKSKDYIN
 ASYIRIVNCGEEYFYIATQGPLLSTIDDFWQMVLENNSNVIAMI TREIEGGI IKCYHYWPI SLKKPLELK
 HFRVFL ENYQILQYFIIRMFQVVEKSTGTSHSV KQLQFTKWPDHGTPASADSF IKYIRYARKSHLTGPMV
 VHCSAGIGRTGVFLCVDVVFCAIVKNC SFNIMDI VAQMREQRSGMVQTKEQYHFCYDIVLEVLRKLLTLD


SGPTRRRLE - GFP Tag - V

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_001042357

ORF Size:	1260 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001042357.4
RefSeq Size:	2965 bp
RefSeq ORF:	1263 bp
Locus ID:	26095
UniProt ID:	Q4JDL3
Cytogenetics:	10q11.22
Protein Families:	Druggable Genome, Phosphatase
Gene Summary:	The product of this gene belongs to the family of classical tyrosine-specific protein tyrosine phosphatases. Many protein tyrosine phosphatases have been shown to regulate fundamental cellular processes. The encoded protein appears to be targeted to sites of actin polymerization. A pseudogene of this gene has been defined on chromosome 10. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2014]