

Product datasheet for **RG215753**

MAGIX (NM_001099681) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MAGIX (NM_001099681) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MAGIX
Synonyms:	JM10; PDZX
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG215753 representing NM_001099681 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGAGCCGCGCACAGGGGGCGCCGCAACCCTAAGGGGAGCAGAGGAAGCCGCGGCCCTCCCGCTCG
CGGGCCCTAGCGCCCGGAGCTCCTGGCGCGTTGGACGCGCGCCCCCTGGCGGCGCGAGCTGCGGTGCA
CGTGGCAGCGCTGGTACGCAGGGCGGGCGCCACATTGCGCCTGCGCCGGAAGGAGGCTGTTAGCGTGCTG
GACTCTGCGGACATAGAGGTACAGACAGTCGCCTGCCTCATGCCACTATTGTGGATCACCAGCCCCAGG
TCGGGGACCTCGTGCTCCACATCAACGGAGAGTCAACGCAGGGCCTCACCCATGCCAGGCCGTGGAGCG
GATCCGAGCTGGAGGCCCCAGCTCCACCTGGTTATTCGTCGGCCTCTGGAGACCCACCCTGGCAAGCCT
CGAGGGGTGGGAGAGCCCCGAAAAGGAGTTGTCCCGTCATGGCCAGATCGCAGCCAGATCCTGGAGGGC
CGGAGGTAACGGGGTCTCGCAGCAGCAGCACTCCCTAGTTTACGACCCTCCATCCCGGACGACGCTCAA
GAAGACCCGGGGCAGCCCGGAGCCTAGTCCAGAGGCGGCCGATGGCCCCACGGTTTCTCCTCTGAG
CGCCGCGCTGAGGATCCCAACGACCAGATCCCGGGTCCCGGGGCCCTGGCTAGTGCCAGCGAGGAAC
GGCTCTCGCGGGCCCTAGGGTCCGGGGGCGAGCGAGCTCGCTCAGGAGATGGCAGCCGGAAGGCGGAG
ACAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >RG215753 representing NM_001099681
 Red=Cloning site Green=Tags(s)

MEPRTGGAANPKGSRGSRGPSPLAGPSARQLLARLDARPLAARAAVDVAALVRRAGATLRLRRKEAVSVL
 DSADIEVTDSRLPHATIVDHRPQVGDVLHINGESTQGLTHAQAVIRAGGPQLHLVIRPLETHPGKP
 RGVGEPRKGVVPSWPDSPDPGGPEVTGSRSSSTSLVQHPPSRTLKKTGSPSPSEAAADGPTVSPPE
 RRAEDPNDQIPGSPGPWLVPSEERLSRALGVRGAAQLAQEMAAGRHH

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001099681

ORF Size: 774 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001099681.1](#), [NP_001093151.1](#)

RefSeq Size: 1922 bp


RefSeq ORF: 777 bp

Locus ID: 79917

UniProt ID: [Q9H6Y5](#)

Cytogenetics: Xp11.23

Product images:



Circular map for RG215753