

Product datasheet for **RG215846**

MAGIX (NM_001099682) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MAGIX (NM_001099682) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MAGIX
Synonyms:	JM10; PDZX
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG215846 representing NM_001099682 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGAGCCGCGCACAGGGGGCGCCGCAACCCTAAGGGGAGCAGAGGAAGCCGCGGCCCTCCCCGCTCG
CGGGCCCTAGCGCCCGGCAGCTCCTGGCGCGTTGGACGCGCGCCCCCTGGCGGCGCAGCTGCGGTGCA
CGTGGCAGCGCTGGTACGCAGGGCGGGCGCCACATTGCGCCTGCGCCGGAAGGAGGCTGTTAGCGTGCTG
GACTCTGCGGACATAGAGGTACAGACAGTCGCCTGCCTCATGCCACTATTGTGGATCACCAGCCCCAGG
TCGGGGACCTCGTGCTCCACATCAACGGAGAGTCAACGCAGGGCCTACCCATGCCAGGCCGTGGAGCG
GATCCGAGCTGGAGGCCCCAGCTCCACCTGTTATTCTGTCGGCCTCTGGAGACCCACCTGGCAAGCCT
CGAGGGGTGGGAGAGCCCCGAAAAGGAGTTGATCGCAGCCAGATCCTGGAGGGCCGGAGGTAACGGGGT
CTCGCAGCAGCAGCACTTCCCTAGTTTCAGCACCTCCATCCCGGACGACGCTCAAGAAGACCCGGGGCAG
CCCGGAGCCTAGTCCAGAGGGCGGCCGCGATGGCCCCACGGTTTCTCCTCCTGAGCGCCGCGCTGAGGAT
CCCAACGACCAGATCCCGGGTTCCCGGGGCCCTGGCTAGTGCCAGCGAGGAACGGCTCTCGCGGGGCC
TAGGGTCCGGGGGCGAGCGAGCTCGCTCAGGAGATGGCAGCCGGAAGGCGGAGACAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >RG215846 representing NM_001099682
 Red=Cloning site Green=Tags(s)

MEPRTGGAANPKGSRGSRGPSPLAGPSARQLLARLDARPLAARAAVDVAALVRRAGATLRLRRKEAVSVL
 DSADIEVTDSRLPHATIDHRPQVGDVLHNGESTQGLTHAQAVIRAGGPQLHLVIRPLETHPGKP
 RGVGEPRKGVDRSPDPGGPEVTGSRSSSTSLVQHPPSRTTLKKTRGSPEPSPEAAADGPTVSPERRAED
 PNDQIPGSPGPWLVPSEERLSRALGVRGAAQLAQEMAAGRHH

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001099682

ORF Size: 759 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001099682.1](#), [NP_001093152.1](#)

RefSeq Size: 1907 bp


RefSeq ORF: 762 bp

Locus ID: 79917

UniProt ID: [Q9H6Y5](#)

Cytogenetics: Xp11.23

Product images:



Circular map for RG215846