

Product datasheet for **RG215871**

Angiotensin Converting Enzyme 1 (ACE) (NM_152830) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Angiotensin Converting Enzyme 1 (ACE) (NM_152830) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Angiotensin Converting Enzyme 1
Synonyms:	ACE1; CD143; DCP; DCP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG215871 representing NM_152830
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGCCAGGGTTGGGCTACTGCAGGACTTCCCAGCCTCCTCTTCTGCTGCTCTGCTACGGGCACCTC
 TGCTGGTCCCAGCCAGGAGGCATCCCAACAGGTGACAGTCACCCATGGGACAAGCAGCCAGGCAACAAC
 CAGCAGCCAGACAACCACCCACCAGGCGACGGCCACCAGACATCAGCCAGAGCCAAAACCTGGTGACT
 GATGAGGCTGAGGCCAGCAAGTTTGTGGAGGAATATGACCGGACATCCCAGGTGGTGTGGAACGAGTATG
 CCGAGGGCAACTGGAACAACACCAACATCACCACAGAGACCAGCAAGATTCTGCTGCAGAAGAACAT
 GCAAATAGCCAACCACACCCTGAAGTACGGCACCCAGGCCAGGAAGTTTGTATGTGAACAGTTGCAGAAC
 ACCACTATCAAGCGGATCATAAAGAAGTTTCAGGACCTAGAACGGGCAGCACTGCCTGCCAGGAGCTGG
 AGGAGTACAACAAGATCCTGTTGGATATGGAACACCTACAGCGTGGCCACTGTGTGCCACCCGAATGG
 CAGCTGCCTGCAGCTCGAGCCAGATCTGACGAATGTGATGGCCAGTCCCGGAAATATGAAGACCTGTTA
 TGGGCATGGGAGGGCTGGCGAGACAAGGCGGGGAGAGCCATCCTCCAGTTTTACCCGAAATACGTGGAAC
 TCATCAACCAGGCTGCCCGGCTCAATGGCTATGTAGATGCAGGGGACTCGTGGAGGTCTATGTACGAGAC
 ACCATCCCTGGAGCAAGACCTGGAGCGGCTCTCCAGGAGCTGCAGCCACTCTACCTCAACCTGCATGCC
 TACGTGCGCCGGGCCCTGCACCGTCACTACGGGGCCAGCACATCAACCTGGAGGGGCCATTCTGCTC
 ACCTGTGGGGAACATGTGGGCGCAGACCTGGTCCAACATCTATGACTTGGTGGTGGCCCTCCCTTCAGC
 CCCCTCGATGGACACCACAGAGGCTATGCTAAAGCAGGGCTGGACGCCAGGAGGATGTTAAGGAGGCT
 GATGATTTCTTACCTCCCTGGGGCTGTGCCCGTGCCTCCTGAGTTCTGGAACAAGTCGATGCTGGAGA
 AGCCAACCGAGCGGGCGGAGGTGGTCTGCCACGCCTCGGCCTGGGACTTCTACAACGGCAAGGACTCCG
 GATCAAGCAGTGCACCACCGTGAACCTGGAGACCTGGTGGTGGCCACCACGAAATGGGCCACATCCAG
 TATTTTCATGCAGTACAAAGACTTACCTGTGGCCTTGAGGGAGGGTGCCAACCCCGCTTCCATGAGGCCA
 TTGGGGACGTGCTAGCCCTCTCAGTGTCTACGCCAACGACCTGCACAGTCTCAACCTGCTGAGCAGTGA
 GGGTGGCAGCGACGAGCATGACATCAACTTTCTGATGAAGATGGCCCTTGACAAGATCGCCTTTATCCCC
 TTCAGCTACCTCGTCGATCAGTGGCGCTGGAGGGTATTTGATGGAAGCATCACCAGGAGAACTATAACC
 AGGAGTGGTGGAGCCTCAGGCTGAAGTACCAGGGCCTTGCCCCCAGTGCCAGGACTCAAGGTGACTT
 TGACCCAGGGGCAAGTCCACATTCCTTCTAGCGTGCCTTACATCAGGTAATTTGTGAGCTTCATCATC
 CAGTTCCAGTTCCACGAGGCACTGTGCCAGGCAGCTGGCCACACGGGCCCTGCACAAGTGTGACATCT
 ACCAGTCCAAGGAGGCCGGCAGCGCCTGGCGACCGCCATGAAGCTGGGCTTCAGTAGGCCGTGGCCGGA
 AGCCATGCAGCTGATCACGGGCCAGCCAAACATGAGCGCCTCGGCCATGTTGAGTACTTCAAGCCGCTG
 CTGGACTGGCTCCGCACGGAGAACGAGCTGCATGGGGAGAAGCTGGGCTGGCCGAGTAACTGGACGC
 CGAACTCCGCTCGCTCAGAAGGGCCCTCCAGACAGCGCCGCGTCAGCTTCTGGGCCTGGACCTGGA
 TGCGCAGCAGGCCCGCTGGGCCAGTGGCTGCTGCTTCTGGGCATCGCCCTGCTGGTAGCCACCCTG
 GGCTCAGCCAGCGGCTTTCAGCATCCGCCACCGCAGCCTCCACCGCACTCCACGGGCCAGTTCCG
 GCTCCGAGGTGGAGCTGAGACTCC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA


Protein Sequence: >RG215871 representing NM_152830
 Red=Cloning site Green=Tags(s)

MGQGWATAGLPSLLFLLLCYGHPLLVPSEQEASQVTVTHGTSSQATTSSQTTTHQATAHQ TSAQSPNLVT
 DEAEASKFVEEYDRTSQVWVNEYAEANWNYNTNITTETSKILLQKNMQIANHTLK YGTQARKFDVNQLQN
 TTIKRIIKKVQDLERAALPAQELEEYNKILLDMETTYSVATVCHPNGSCLQLEPDLTNVMATSRKYEDLL
 WAWEGWRDKAGRAILQFYPKYVELINQAARLNGYVDAGDSWRSMYETPSLEQDLERLFQELQPLYLNLHA
 YVRRALHRHYGAQHINLEGPPIPAHLLGNMWAQTWSNIYDLVVPFAPSMDTTEAMLKQGWTPRRMFKEA
 DDFFTSLGLLPVPPEFWNKSMLEKPTDGREVVCHASAWDFYNGKDFRIKQCTTVNLEDLVVAHHEMGIQ
 YFMQYKDLPVALREGANPGFHEAIGDVLALSVSTPKHLHSLNLLSSEGGSD EHDINFLMKMALDKIAFIP
 FSYLVDQWRWRVFDGSITKENYNQEWLRLKYQGLCPPVPRTQGDFDPGAKFHIPSSVPIRYFVSFII
 QFQFHEALCQAAGHTGPLHKCDIYQSKEAGQLATAMKLGFSRPWPEAMQLITGQP NMSASAMLSYFKPL
 LDWLRTENELHGEKLGWPQYNWTPNSARSEGPLPDSGRVSFLGLDLDAQQARVQWLLFLGIALLVATL
 GLSQRLF SIRHRS LHRHSHGPQFGSEVELRHS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_152830

ORF Size: 2196 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_152830.3](#)

RefSeq Size: 3264 bp

RefSeq ORF: 2199 bp

Locus ID: 1636

UniProt ID: [P12821](#)


Cytogenetics: 17q23.3

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Protease, Secreted Protein, Transmembrane

Protein Pathways: Hypertrophic cardiomyopathy (HCM), Renin-angiotensin system

Gene Summary: This gene encodes an enzyme involved in blood pressure regulation and electrolyte balance. It catalyzes the conversion of angiotensin I into a physiologically active peptide angiotensin II. Angiotensin II is a potent vasopressor and aldosterone-stimulating peptide that controls blood pressure and fluid-electrolyte balance. This angiotensin converting enzyme (ACE) also inactivates the vasodilator protein, bradykinin. Accordingly, the encoded enzyme increases blood pressure and is a drug target of ACE inhibitors, which are often prescribed to reduce blood pressure. This enzyme additionally plays a role in fertility through its ability to cleave and release GPI-anchored membrane proteins in spermatozoa. Many studies have associated the presence or absence of a 287 bp Alu repeat element in this gene with the levels of circulating enzyme. This polymorphism, as well as mutations in this gene, have been implicated in a wide variety of diseases including cardiovascular pathophysiologies, psoriasis, renal disease, stroke, and Alzheimer's disease. Regulation of the homologous ACE2 gene may be involved in progression of disease caused by several human coronaviruses, including SARS-CoV and SARS-CoV-2. Alternative splicing results in multiple transcript variants encoding both somatic (sACE) and male-specific testicular (tACE) isoforms. [provided by RefSeq, Sep 2020]

Product images:



Circular map for RG215871