

Product datasheet for **RG215907**

UEVLD (NM_018314) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	UEVLD (NM_018314) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	UEVLD
Synonyms:	ATTP; UEV3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG215907 representing NM_018314 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAGTTCGACTGCGAGGGCCTGAGACGGCTGCTTGGCAAGTACAAGTTCAGGGACCTAACTGTGGAAG
 AACTAAGGAATGTAAATGTATTTTCCACATTTCAAATATTCCATGGACACCTATGTTTTAAAGATAG
 TTCTCAGAAAGACCTGCTGAATTTTACTGGCACAATTCCTGTGATGTATCAGGGTAATACATATAACATA
 CCAATTCGTTTCTGGATTTTGGATTCTCACCTTTTCGCTCCCCCTATTTGCTTCTTGAAGCCAAGTCAA
 ATATGGGAATCTTAGTCGAAAACATGTGGATGCTCAAGGCAGAAATATATTTGCCCTATCTCCAAAAGT
 GAGCCATCCTAAATCTGTCAATTGTTGGATTAATTAAGAAATGATTGCCAAGTTTCAAGAGGAAGTCCC
 ATGTATTTCTATCATCATCTGATGAGGCACGGCAGGTAGACTTGCTAGCCTATATTGCAAAAATCACTG
 AAGGTGTTTCAGATACAAATTCAGAGCTGGGCAATCATGAGAATAAAACAGTCAATAAAATTAAGTGT
 GGTTGGAGGTGGAGAAGTGGTATTGCCTGCACATTAGCAATTTAGCAAAAGGTATTGCAGACAGGCTT
 GTCCTCTTAGACCTCTCAGAAGGGACTAAAGGAGCCAGATGGACCTTGAATCTTCAACCTTCCTAATG
 TGGAGATCAGCAAAGATTTGCTGCCTCTGCTCATTCCAAGGTGGTGATCTTACAGTCAACTCTTTGGG
 TAGTTCTCAGTCGTACCTTGATGTGGTACAGAGCAATGTGGATATGTTTACAGACCTTGTCCAGCTCTG
 GGACATTATAGTCAACACAGTGTCTGCTGCTGATCTCAACCAGTGGAAATCATGACCTATGTAACAT
 GGAAACTGAGTACATTTCTGCAAATCGAGTGATCGGAATTGGATGTAATCTGGATTACAGAGATTACA
 GTATATTATTACAAATGTTTGAAGGCACAGACTTCAGGCAAAGAAGTATGGGTTATTGGCGAGCAAGGA
 GAAGACAAAGTGCTCACATGGAGTGGCCAAAGAAGTAGTGAGTCTTACCTCTCAAGTGCAGCTGTCCA
 ACAGGGATATTATGATA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


[View online »](#)

Protein Sequence: >RG215907 representing NM_018314
 Red=Cloning site Green=Tags(s)


MEFDCEGLRRLGKYKFRDLTVEELRNVNFFPHFKYSMDTYVFKDSSQKDLLNFTGTIPVMYQGNTYNI
 PIRFWILDSPFAPPICFLKPTANMGILVGKHVDAQGRIYLPYLQNWHPKSVIVGLIKEMIAKFQEELP
 MYSLSSSDEARQVDLLAYIAKITEGVSDTNSKSWANHENKTVNKITVVGGEGLIACTLAISAKGIADRL
 VLLDLSEGTKGATMDLEIFNLPNVEISKDLSASAHSKVVIPTVNSLGSSQSYLDVVQSNVDMFRALVPAL
 GHYSQHSVLLVASQPVEIMTYVTWKLSTFPANRVIGIGCNLDSQRLQYIITNVLKAQTSKKEVWVIGEQQ
 EDKVLTSWGQEEVVSLSQVQLSNRDI

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_018314

ORF Size: 1137 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_018314.2](#), [NP_060784.2](#)

RefSeq Size: 1984 bp

RefSeq ORF: 1140 bp

Locus ID: 55293


UniProt ID: [Q8IX04](#)

Cytogenetics: 11p15.1

Domains: Idh

Gene Summary: Possible negative regulator of polyubiquitination.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG215907