

Product datasheet for **RG215925**

TSGA10IP (NM_152762) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	TSGA10IP
Synonyms:	FAM161C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



ORF Nucleotide Sequence: >RG215925 representing NM_152762
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGGGCAGGACACCGATATGCTAAATACCTACCAACAGTTGGTTAGGACCCCGTCGGTGCGACACAGGGC
 AGGACGTGCGGCTCCAGGCTCCAGGAACAGAACGGGGCTGCTCAAGCTGCTGTCGACCGTCTCTCAGGA
 CAAGCAGGGCTGCCTGGGAGTGGTGATGGTGTCCAAATCAAGACCTGCAGCAGAGGTCTCAGAGCTCA
 AGGCAGACAGCAAGAAGGACCGCAAGCCAGGGGCCAGAGCAAGAAAGACAGGGCTCTGAAGAGTCTG
 AAGATCACTTCCCGCTTCTTCTCGGAAACCTCCTTCCCTTCCAGTGGGCTGGGAGAGCATCGCCAC
 AGATGTCCGGGTGTGCTTCAGCCAAGCTCCCAACCCCTGGCCACCAAGCCCTGCCATGCCCTCCTCG
 TTCTCCAGCGTCAGTCCAGGCGAAGTCCACGGCCAACCTCCAGAGGCCATGGTGCTGCTGGAAGA
 CAGAGGCGCAAAACCTGAAGGCGAGACAGCAGCTGGGAGCCTGGGCGGTGTCTCCATCCCTACTGGCAA
 AGGGGAGCTAGGATCAGAGCCCCCAGTGGTCTCCAGCTGCCTGGGAGGAGGCCAGGATCAGGATCGGCG
 TCCGACAAGCAGGTCCAGCTGCAAAGCCTGGGTGCTGAGGAGGCCGAGAGGGGTCTGAGTTCTGGGGTGC
 TGCCCCAGCGCCCCAGGAGGGGTGCTATCTCAGAGGAGGAGCAATTTTCAGAGGCCACAGAGGAGGTGA
 GGAGGGAGAGCACAGGACTCCCTGCAGAAGGAGGGTGGTTGTGTCAGAGAAAGGGGAGATTTCTGGAGAG
 GAAGCCTCCGATGAGGGAGAAGTGCAGGGCCAGAGCCAGGGGAGCAGCCCCAGCTTCAACAACCTCCGAA
 GGCGACAATGGAGGAAGACAAGGGCCAAGGAGCTGCAGGGGCCATGGGACCTGGAGAAGCTGCACAGGCA
 GCTACAGAGAGACCTGGACTGTGGCCCCAAAAGCTGCCCTGGAAGACTTTGAGGGCTGCCTTCCAGGCC
 TCCAAGCGGAATGGAAGGCCATGCTCGGGATACGATGAACTTTCTGTCTGCCAACCTCCCTAATC
 GCACCTTCCACAACGACAGGAAGCCACCAGGAGCCTGCTGCAGGCCTGGGAGCGGCAGCGCAGGAGGA
 GCGGCAGCAGGCCGAGCTGCGGCGGGCCCGACACAGCATGTACAGCGGCAGGTGGCCCACTGCCTGGCA
 GCCTACGCACCCAGAGGAGCCGGGGCCCTGGGGCGGCCAGCGCAAGCTGGAGGAGCTGAGGCGCCAGG
 AGCGACAGCGCTTTGCTGAGTACCAGGCGGAGCTGCAAGGCATCCAGCAGAGGTGCAGGCCCGGCCCTT
 CCTGTTCCAGCAGGCTATGCAGGCCAATGCCGGCTCACCCTCACTCGGCGCTTCTCCAGGTGCTGTCA
 GCACTGGGGCTGGATGAGGAGCAGCTGCTGTCTGAGGCAGGAAAGGTGGACAGAGAGGGCACCCCCAGGA
 AACCCAGGAGCCACAGGTCACTGGGGGTGAGAATGGAGCACTCTCCTCAGAGGCCCCCAAGGACAGAACC
 CACCGGCAGCCAGCCTGACAGGCACTACAACCCAGCCTGGACCCGGAGTGCAGTCCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG215925 representing NM_152762
 Red=Cloning site Green=Tags(s)

MGQDTMLNTYQQLVVRTPSVRPGQDVRLQAPGTRTGLLKLSTVSQDKQGCLGSGDGVPNQDLQQRSQSS
 RQTAKKDRKPRGQSKKGQSESEDFPLLPKPSFPFQWAWESIATDVRAVLQPSSTPGHQALPMPSS
 FSQRQSRRKSTANLPEAHGCCWKTEAQNLIKARQLGAWGGVSIPTGKGLGSEPPSGLQLPGRRPGSGSA
 SDKQVQLQSLGAEEAERGLSSGVLPRPRRGSISEEEQFSEATEEAEEGEHRTPCRRRAGCQRKGQISGE
 EASDEGEVQGSQSSPSFNLRQRWRKTRAKELQGPWDLEKLHRQLQRDLDCGPQKLPWKTLLRAAFQA
 SKRNGKAYASGYDETFVSANLPNRTFHKRQEATRSLLQAWERQRQEERQQAELRRARTQHVQRQVAHCLA
 AYAPRGRGPGAAQRKLEELRRQERQRF AEYQAE LQGIQHRVQARPFLFQQAMQANARLTVTRRF SQVLS
 ALGLDEEQLLSEAGKVDREGTPRKPRSHRSVGRMEHSPQRPPRTEPTGSQPD RHYNPSLDPECSF

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_152762

ORF Size: 1668 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_152762.1](#), [NP_689975.1](#)

RefSeq Size: 1726 bp

RefSeq ORF: 1671 bp
 Locus ID: 254187
 UniProt ID: [Q3SY00](#)
 Cytogenetics: 11q13.1

Product images:



Circular map for RG215925