

Product datasheet for **RG215940**

14-3-3 beta (YWHAB) (NM_003404) Human Tagged ORF Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | 14-3-3 beta (YWHAB) (NM_003404) Human Tagged ORF Clone |
| Tag: | TurboGFP |
| Symbol: | 14-3-3 beta |
| Synonyms: | GW128; HEL-S-1; HS1; KCIP-1; YWHAA |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-AC-GFP (PS100010) |
| E. coli Selection: | Ampicillin (100 ug/mL) |
| ORF Nucleotide Sequence: | >RG215940 representing NM_003404 Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACAATGGATAAAAGTGAGCTGGTACAGAAAGCCAACTCGCTGAGCAGGCTGAGCGATATGATGATA
TGGCTGCAGCCATGAAGGCAGTCACAGAACAGGGGCATGAACTCTCCAACGAAGAGAGAAATCTGCTCTC
TGTTGCCATAAGAATGTGGTAGGCGCCCGCGCTCTTCTGGCGTGTCTCCAGCATTGAGCAGAAA
ACAGAGAGGAATGAGAAGAAGCAGCAGATGGGCAAAGAGTACCGTGAGAAGATAGAGGCAGAAGTGCAGG
ACATCTGCAATGATGTTCTGGAGCTGTTGGACAAATATCTTATTCCAATGCTACACAACCAGAAAGTAA
GGTGTCTACTTGAAAATGAAAGGAGATTATTTTAGGTATCTTTCTGAAGTGGCATCTGGAGACAACAAA
CAAACCACTGTGTCGAACTCCCAGCAGGCTTACCAGGAAGCATTGAAATTAGTAAGAAAGAAATGCAGC
CTACACACCCAATTCGTCTTGGTCTGGCACTAAATTTCTCAGTCTTTACTATGAGATTCTAAACTCTCC
TGAAAAGGCTGTAGCCTGGCAAAAACGGCATTGATGAAGCAATTGCTGAATTGGATACGCTGAATGAA
GAGTCTTATAAAGACAGCACTCTGATCATGCAGTTACTTAGGGACAATCTCACTCTGTGGACATCGGAAA
ACCAGGGAGACGAAGGAGACGCTGGGGAGGGAGAGAAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG215940 representing NM_003404
 Red=Cloning site Green=Tags(s)

MTMDKSELVQKAKLAEQAERYDDMAAAMKAVTEQGHLSNEERNLLSVAYKNVVGARRSSWRVISSIEQK
 TERNEKKQQMGKEYREKIEAELQDICNDVLELLDKYLIPNATQPESKVFYLMKMGDYFRYLSEVASGDNK
 QTTVSNSQQAYQEAFAEISKEMQPTHPIRLGLALNFSVFYYEILNSPEKACSLAKTAFDEAIAELDTLNE
 ESYKDSTLIMQLLRDNLTLWTSENQGDGEGDAGEGEN

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_003404

ORF Size: 738 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_003404.5](#)

RefSeq Size: 3110 bp

RefSeq ORF: 741 bp

Locus ID: 7529

UniProt ID: [P31946](#)

Cytogenetics: 20q13.12


Domains: 14-3-3

Protein Families: Druggable Genome

Protein Pathways: Cell cycle, Neurotrophin signaling pathway, Oocyte meiosis

Gene Summary: This gene encodes a protein belonging to the 14-3-3 family of proteins, members of which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals. The encoded protein has been shown to interact with RAF1 and CDC25 phosphatases, suggesting that it may play a role in linking mitogenic signaling and the cell cycle machinery. Two transcript variants, which encode the same protein, have been identified for this gene. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG215940