

## Product datasheet for **RG216009**

### **FAM177A1 (NM\_173607) Human Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** FAM177A1 (NM\_173607) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** FAM177A1  
**Synonyms:** C14orf24  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG216009 representing NM\_173607  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGAAGTGGGCTTACCGGCCATTACCCTCTTTCTCACCAGCGCCAGCAGCCCTGTGGTGGCGACGACGA  
TGGACCAGGAGCCAGTGGGCGGTGTGGAACGAGGAGAAGCCGTCGCAGCCTCGGGAGCTGCGGCCCGCC  
GGCATTGGGGAATCTGCAGGGCAGATGAGTAACGAAAGAGGCTTTGAAAATGTAGAACTGGGAGTCATA  
GGAAAAAAGAAGAAAGTCCCAAGGAGAGTCATCCACTTTGTTAGTGGTGAACAATGGAAGAATATAGCA  
CAGATGAAGACGAAGTTGATGGCCTGGAGAAGAAAGATGTTTTGCCTACTGTTGATCCGACAAAATTAC  
CTGGGGTCCCTACTTATGGTTTTACATGCTTCGGGCTGCTACATCAACTCTCAGTGTGTGACTTCCTT  
GGAGAGAAGATTGCATCTGTTTTGGGTATCAGCACCCCAAAGTACCAATATGCCATTGATGAATATTATC  
GGATGAAGAAGGAGGAAGAAGAAGAAGAAGAAAACAGGATGTCTGAAGAAGCAGAAAAACAATATCA  
ACAGAATAAATTGCAGACTGATTCCATTGTTTCAGACAGATCAACCAGAGACAGTGATATCCAGCTCATT  
GTGAATGTCAATTTTGAATGGAGGGAGACAGTGAAGTAATTATGGAAAGCAAGCAAAATCCAGTCTCTG  
TCCCACCA

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG216009 representing NM\_173607  
Red=Cloning site Green=Tags(s)

MEVGLPAITLFLTSASSPVVATTMDQEPVGGVERGEAVAASGAAAAAFGESAGQMSNERGFENVELGVI  
 GKKKKVPRRVIHFVSGETMEEYSTDEDEVDGLEKDDVLPVDPKLTWGPYLFYMLRAATSTLSVCDLF  
 GEKIASVLGISTPKYQYAIDEYRMMKEEEEEEEEEENRMSEEAKEQYQNKLQTDTSIVQTDQPETVISSSF  
 VNVNFEMEGDSEVIMESKQNPVSVPP

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**Plasmid Map:**



**ACCN:** NM\_173607

**ORF Size:** 708 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_173607.5</a>
<b>RefSeq Size:</b>	2946 bp
<b>RefSeq ORF:</b>	711 bp
<b>Locus ID:</b>	283635
<b>UniProt ID:</b>	<a href="#">Q8N128</a>
<b>Cytogenetics:</b>	14q13.2
<b>Gene Summary:</b>	This gene encodes a member of a conserved protein family. Alternative splicing results in multiple transcript variants. This gene is thought to be associated with susceptibility to juvenile idiopathic arthritis. [provided by RefSeq, Apr 2017]