

Product datasheet for **RG216125**

TM2D2 (NM_001024381) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	TM2D2
Synonyms:	BLP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG216125 representing NM_001024381
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGCACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGAAGAAATGTCTTTTGGCCGTTTGTATTACGTGCATGCAAACAGCGATTTGCAAAGACCGTATGATGA
TGATCATGATCTTACTGGTGAATTACAGACCTGATGAATTTATAGAATGTGAAGACCCAGTGGATCATGT
TGGAATGCAATGCATCCCAGGAACCTGGTTATGGTTGTCTCAAGTTCGGCGGTGAGGCTACAGCGAC
GTGGAACACACTTCAGTCCAGTGCCATGCCTAGATGGAATTGAGTGTGCCAGTCTAGGACCTTTCTAC
GAGAAAATAAACCTTGTATAAAGTATACCGGACACTACTTCATAACCACTTTACTCTACTCCTTCTCTCT
GGGATGTTTTGGTGTGGATCGATTCTGTTTGGGACACACTGGCACTGCAGTAGGGAAGCTGTTGACGCTT
GGAGGACTTGGGATTTGGTGGTTTGTGACCTATTTTGCTAATTACTGGAGGGCTGATGCCAAGTGATG
GCAGCAACTGGTGCAGTGTTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG216125 representing NM_001024381
Red=Cloning site Green=Tags(s)

MKKCLLPVLITCMQTAICKDRMMIMILLVNYRPDEFIECEDPVDHVGNGATASQELGYGCLKFGGQAYS
VEHTSVQCHALDIECASPRFLRENKPKIKYTGHYFITLLYSFFLGCFGVDRFLGHTGTAVGKLLTL
GGLGIWWFVDLILLITGGLMPDGSNWCTVY

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: NM_001024381

ORF Size: 513 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001024381.2](#)

RefSeq Size:	3334 bp
RefSeq ORF:	516 bp
Locus ID:	83877
UniProt ID:	Q9BX73
Cytogenetics:	8p11.22
Protein Families:	Druggable Genome, Transmembrane

Gene Summary: The protein encoded by this gene contains a structural module related to that of the seven transmembrane domain G protein-coupled receptor superfamily. This protein has sequence and structural similarities to the beta-amyloid binding protein (BBP), but, unlike BBP, it does not regulate a response to beta-amyloid peptide. This protein may have regulatory roles in cell death or proliferation signal cascades. This gene has multiple alternatively spliced transcript variants which encode two different isoforms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG216125