

Product datasheet for **RG216133**

PMS2L5 (NM_174930) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	PMS2L5
Synonyms:	MGC34222; PMS7; postmeiotic segregation increased 2-like 5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG216133 representing NM_174930
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTGGGGTAGAAGAAGAAACTTCGAAGGCTTAATGATGTCACCATTTCTACCTGCCACGTATCGGCGA
AGGTTGGGACTCGACTGGTGTGGATCACGATGGGAAATCATCCAGAAACCCCTACCCCAACCCAG
AGGGACCACAGTCAGCGTGAAGCAGTTATTTCTACGCTACCTGTGCGCCATAAGGAATTTCAAAGGAAT
ATTAAGAAGAAACGTGCCTGCTTCCCTTCGCCTTCTGCCGTGATTGTCAGTTCCTTGAGGGCTCCCGAG
CCATGCTTCTGTACAGCCTGCAAACTGACTCCTAGAGTACCCCAACCCCTGCTCCTTGAGGGA
CAACGTGATCACTGTATTGAGCTCTGTCAAGAATGGTCCAGGTTCTTCTAGA

AG**CGGACCG**ACGCGTACGCGCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG216133 representing NM_174930
Red=Cloning site Green=Tags(s)


MWGRRRKLRLNDVTISTCHVSAKVGTRLVFDHDKIIQKTPYPHPRGTTVSVKQLFSTLPVRHKEFQRN
IKKKRACFPFAFCRDCQFLEGSPAMLPVQPAKLTPRSTPPHPCSLIEDNVITVFSSVKNPGSSR

SGPTRRRLE - GFP Tag - V

Restriction Sites: SgfI-RsrII



Cloning Scheme:



ACCN: NM_174930

ORF Size: 402 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_174930.3

RefSeq Size: 765 bp


RefSeq ORF: 404 bp

Locus ID: 5383

Cytogenetics: 7q11.23

Protein Families: Druggable Genome

Product images:



Circular map for RG216133