

## Product datasheet for **RG216172**

### **C2orf80 (NM\_001099334) Human Tagged ORF Clone**

#### Product data:

|                           |                         |
|---------------------------|-------------------------|
| <b>Product Type:</b>      | Expression Plasmids     |
| <b>Tag:</b>               | TurboGFP                |
| <b>Symbol:</b>            | C2orf80                 |
| <b>Synonyms:</b>          | GONDA1                  |
| <b>Mammalian Cell</b>     | Neomycin                |
| <b>Selection:</b>         |                         |
| <b>Vector:</b>            | pCMV6-AC-GFP (PS100010) |
| <b>E. coli Selection:</b> | Ampicillin (100 ug/mL)  |

**ORF Nucleotide Sequence:** >RG216172 representing NM\_001099334  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGAAGAAGGCTCATAAAGAAGGAAATGAAAAAGCTCTTGGGAGATTATATTGGCATCAGACTTCGGG  
AAAATGAATTTGACCCAAAAGGACGGCAACTCACCTTTCTAGATGATATGGCACACTATGACTTGCC  
CATCAGTGTGCTTTGCAATGGCTGGATCCCTCAGAAGACTTAACCTGGCTGGAGTGGGAGGAAGTAAA  
ATACCACTCCATGGCAGACCCATATATCCAATCGTAGAGAACGAGAAGCTATGATTTTATCATCTTATG  
CTGGAATCTTAATGAACAGTATCCCGATTGAGGAAGTCTTTAAAAATTTATGGGGCTGATTCTTCTGCCGA  
TTCTGGTACCATCAAGGTTCCCGAGTTTCATCTCTCCGCCTCTCCTTGACCCCTTTGCCATGTTAACA  
GCACCCAAAGCAGCAGCATACGCCCGCAACAGAGTGTCAAGTCAAGAAAGGTAACAACAAATAAAATG  
CCACCAGCATCTGCAAAGGAAGCAAATGCCACAGAAATGGAATCATCACAAGGTTTTCAGACACACA  
GCCAAAGCACAAGTCACT

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >RG216172 representing NM\_001099334  
Red=Cloning site Green=Tags(s)


MERRLIKEMKLLGDYIGIRLRENEFDPKGRRQLTFLDDMAHYDLAISVALQWLDPSEDLTWLEWEEK  
IPLHGRPIYPNRREREAMILSSYAGILMNSIPIEEVFKIYGADSSADSGTIKVPRVSSLRLSLHPFAMLT  
APKAAAYARKQSVKSRKVTTNKNATSISAKEANATEWKSSQRFSDTQPKHKVT

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** SgfI-MluI



## Cloning Scheme:



**ACCN:** NM\_001099334

**ORF Size:** 579 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

**RefSeq:** [NM\\_001099334.2](#), [NP\\_001092804.1](#)


**RefSeq Size:** 1229 bp

RefSeq ORF: 582 bp

Locus ID: 389073

UniProt ID: [Q0P641](#)

Cytogenetics: 2q33.3

**Product images:**


Circular map for RG216172