

Product datasheet for **RG216177**

Profilin 2 (PFN2) (NM_002628) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Profilin 2 (PFN2) (NM_002628) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: Profilin 2
Synonyms: D3S1319E; PFL
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG216177 representing NM_002628
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCGGTTGGCAGAGCTACGTGGATAACCTGATGTGCGATGGCTGCTGCCAGGAGGCCGCCATTGTCTG
GCTACTGCGACGCCAAATACGCTCTGGGCAGCCACGGCCGGGGCGCTTTTCAGAGCATTACGCCAATAGA
AATAGATATGATTGTAGGAAAAGACCGGGAAGTTTCTTTACCAACGTTTGACTCTTGGCGGAAGAAA
TGCTCAGTGATCAGAGATAGTCTATACGTCGATGGTACTGCACAATGGACATCCGGACAAAGAGTCAAG
GTGGGGAGCCAACATAACAATGTGGCTGTGGCAGAGCTGGTAGAGCATTGTTATAGTCATGGGAAAGGA
AGGTGTCCACGGAGGCACACTTAACAAGAAAGCATATGAACTCGCTTTATACCTGAGGAGGTCTGATGTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG216177 representing NM_002628
Red=Cloning site **Green**=Tags(s)
MAGWQSYVDNLMCDGCCQEAAIVGYCDAKYVWAATAGGVFQSITPIEIDMIVGKDREGFFTNGLTLGAKK
CSVIRDSLIVDGDCTMDIRTKSQGGEPTYNVAVGRAGRALVIVMKEGVHGGTLNKKAYELALYLRRSDV

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_002628

ORF Size: 420 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_002628.5](#)

RefSeq Size: 1950 bp

RefSeq ORF: 423 bp

Locus ID: 5217
UniProt ID: [P35080](#)
Cytogenetics: 3q25.1
Domains: PROF
Protein Families: Druggable Genome
Protein Pathways: Regulation of actin cytoskeleton
Gene Summary: The protein encoded by this gene is a ubiquitous actin monomer-binding protein belonging to the profilin family. It is thought to regulate actin polymerization in response to extracellular signals. There are two alternatively spliced transcript variants encoding different isoforms described for this gene. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG216177