

## Product datasheet for **RG216267**

### **FGF5 (NM\_033143) Human Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** FGF5 (NM\_033143) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** FGF5  
**Synonyms:** HBGF-5; Smag-82; TCMGLY  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG216267 representing NM\_033143  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGAGCTTGTCTCCTCCTCCTCCTCCTCCTTTCAGCCACCTGATCCTCAGCGCCTGGGCTCACGGGGAGA  
AGCGTCTCGCCCCAAAGGGCAACCCGGACCCGCTGCCACTGATAGGAACCCTAGAGGCTCCAGCAGCAG  
ACAGAGCAGCAGTAGCGCTATGTCTTCTCCTTCTGCCTCCTCCTCCCCGCAGCTTCTCTGGGCAGCCAA  
GGAAGTGGCTTGGAGCAGAGCAGTTTCCAGTGGAGCCCCCTGGGGCGCCGACCGGCAGCCTCTACTGCA  
GAGTGGGCATCGTTTTCCATCTGCAGATCTACCCGGATGGCAAAGTCAATGGATCCACGAAGCCAATAT  
GTTAAGCCAAGTTCACAGA

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >RG216267 representing NM\_033143  
**Red**=Cloning site **Green**=Tags(s)

MSLSFLLLLFFSHLILSAWAHGEKRLAPKGQPGAATDRNPRGSSSRQSSSSAMSSSSASSSPAASLGSQ  
GSGLEQSSSQWSPSGRRTGSLYCRVIGIFHLQIYPDGVKNGSHEANMLSQVHR

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_033143

**ORF Size:** 369 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


**RefSeq:** [NM\\_033143.2](#), [NP\\_149134.1](#)

**RefSeq Size:** 5295 bp

**RefSeq ORF:** 372 bp

**Locus ID:** 2250  
**UniProt ID:** [P12034](#)  
**Cytogenetics:** 4q21.21  
**Protein Families:** Druggable Genome, Secreted Protein  
**Protein Pathways:** MAPK signaling pathway, Melanoma, Pathways in cancer, Regulation of actin cytoskeleton  
**Gene Summary:** The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. This gene was identified as an oncogene, which confers transforming potential when transfected into mammalian cells. Targeted disruption of the homolog of this gene in mouse resulted in the phenotype of abnormally long hair, which suggested a function as an inhibitor of hair elongation. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]

### Product images:



Circular map for RG216267