

Product datasheet for **RG216290**

LCEIB (NM_178349) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Tag: TurboGFP
Symbol: LCEIB
Synonyms: LEP2; SPRL2A
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG216290 representing NM_178349
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCTGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCCTGCCAGCAGAACCAGCAGCAGTGCCAGCCCCCTCCCAAGTGCATCCCCAAGTGCCCTCCCAAGT
GCCTACCCCTAGATGCCCCCAAGTGTCCCCCTAAGTGTCTCCAGTCTCTTCTGCTGCAGTGTCTAG
CTCCGGAGGCTGCTGTGGCTCCAGCTCTGGGGGAAGCTGTGGCTCCAGCTCTGGGGTTGCTGCAGTTCT
GGGGGAGGTGGCTGCTGCCTGAGCCACCACAGGCGCCGTAGGTCCCACTGCCACAGACCCAGAGCTCTG
GCTGCTGCAGCCAGCCCTCCGGGGGCTCCAGCTGCTGTGGAGGAGGAGTGCCAGCACTCTGGAGGCTG
CTGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG216290 representing NM_178349
Red=Cloning site Green=Tags(s)

MSCQQNQQCQPPPKCIPKCPKCLTPRCPPKCPPVSSCCSVSSGGCCGSSSGGSCGSSSGGCCSS
GGGGCCLSHRRRRSHCHRPQSSGCCSQPSGGSSCCGGSGQHSGGCC

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: NM_178349

ORF Size: 354 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_178349.2

RefSeq Size: 1197 bp

RefSeq ORF: 357 bp


Locus ID: 353132

UniProt ID: Q5T7P3

Cytogenetics: 1q21.3

Gene Summary: Precursors of the cornified envelope of the stratum corneum.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG216290