

Product datasheet for **RG216386**

H2AC18 (NM_003516) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: H2AC18 (NM_003516) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: H2AC18
Synonyms: H2A; H2a-615; H2A.2; H2A/O; H2A/q; H2AC19; H2AFO; HIST2H2AA; HIST2H2AA3
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG216386 representing NM_003516
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGGTCTGGCAAGCAAGGAGGCAAGGCCCGCCCAAGCCAAGTCGCGCTCGTCCCAGCTGGCC
TTCAGTCCCAGTAGGGCGAGTGCATCGCTTGTGCGCAAAGGCAACTACGCGGAGCGAGTGGGGCCGG
CGCGCCCGTCTACATGGCTGCGGTCCTCGAGTATCTGACCGCCGAGATCCTGGAGCTGGCGGCAACGCG
GCTCGGGACAACAAGAAGACGCGCATCATCCCTCGTCACCTCCAGCTGGCCATCCGCAACGACGAGGAAC
TGAACAAGCTGCTGGCAAAGTACCATCGCCAGGGCGGCGTCTGCCTAACATCCAGGCCGACTGCT
CCCTAAGAAGACGGAGATCACCAAGGCAAAGGCAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG216386 representing NM_003516
Red=Cloning site Green=Tags(s)
MSGRGKQGGKARAKAKSRSSRAGLQFPVGRVHRLLRKGNYAERVGAGAPVYMAAVLEYLTAEILELAGNA
ARDNKKTRIIPRHLQLAIRNDEELNKLKGVYIAQGGVLPNIQAVLLPKKTESHKAKGK

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_003516.3
RefSeq Size:	534 bp
RefSeq ORF:	393 bp
Locus ID:	8337
UniProt ID:	Q6FI13
Cytogenetics:	1q21.2
Domains:	H2A, histone
Protein Pathways:	Systemic lupus erythematosus
Gene Summary:	Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H2A family. Transcripts from this gene lack polyA tails but instead contain a palindromic termination element. This gene is found in a histone cluster on chromosome 1. This gene is one of four histone genes in the cluster that are duplicated; this record represents the centromeric copy. [provided by RefSeq, Aug 2015]