

Product datasheet for **RG216440**

ZNF557 (NM_024341) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	ZNF557
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PSI00010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG216440 representing NM_024341
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCGGCTGTCGTCTGCCCACTGCCGCTCTGTCTCCCTGTTCCAGCCTCTCAGCGAGAAGGAC
ACACAGAGGGCGGAGAGCTGGTTAATGAGCTCTGAAAAGCTGGCTAAAGGGCTTGGTGACCTTTGAGGA
TGTGGCCGTGGAGTTCACCCAGGAGGAGTGGGCATTGCTGGACCTGCCCAAAGGACACTGTACAGGGAC
GTGATGCTGGAGAAGTGCAGGAACCTGGCCTCACTGGGGAACCAAGTTGATAAACCTAGGCTGATCTCCC
AGCTGGAGCAAGAAGATAAAGTGATGACAGAAGAGAGGAATTCTCTCAGGTACCTGTCCAGATGTGGA
GAATCCATTTAAAGCCAAAGGGTTAACTCCTAAGCTGCATGTTTTTCGAAAAGAACAATCTAGAAATATG
AAAATGGAGAGGAATCATCTTGGAGCAACACTCAACGAATGTAATCAGTGTTTTAAAGTCTTCAGCACAA
AATCTTCCCTTACACGGCACAGGAAGATTCATACTGGAGAAAGACCCTATGGCTGCAGTGAATGTGGAA
ATCCTACAGCAGTAGATCTTACCTTGCTGTTCAAGAGAATCCACAATGGGGAGAAACCTATGAATGC
AATGACTGTGGGAAACCTTCAGCAGCAGATCTTACCTTACTGTTCAAGAGAATCCACAATGGGGAGA
AACCTTACGAATGCAGTGACTGTGGGAAACCTTCAGCAATTCCTCATACCTCAGACCGCACTTGAGAAT
TCACACTGGAGAAAAACGTACAAATGTAACCAAGTGTTCGTGAGTTCGCACTCAGTCAATCTTCACA
AGGCACAAGAGAGTTTATACGGGGGAGGGTCATTATGTATGTAATCAGTGTGGAAGGCTTTCGGCACGA
GGTCATCTCTTTCTCGCACTATAGCATTATACAGGGGAGTACCTTACGAATGCCACGATTGTGGGAG
AACCTTCAGGAGGAGTCTGAATCTGACACAGCACATAAGAAGTCACTGGAGAAAAACCTACACATGT
AATGAGTGTGGAAATCCTTTACCAATAGCTTTTCTCTTACAATTCACAGGAGAATACATAATGGAGAGA
AATCCTATGAGTGCAGTGATTGTGAAAAATCCTTAATGTTCTCTCATCCGTTAAGAAACACATGAGAAC
TCACACTGGAAAAAACCTATGAATGTAATTATGCGGGAAATCCTTCACAAGTAACCTACCTTTCT
GTGCATACGAGAATGCATAATAGGCAATG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG216440 representing NM_024341
Red=Cloning site Green=Tags(s)

MAAVLPPTAALSSLFPASQREGHTEGGELVNELLKSWLKGLVTFEDVAVEFTQEEWALLDPAQRTLYRD
VMLENCRNLASLGNQVDKPRILISQLEQEDKVMTEERGILSGTCPDENPFKAKGLTPKLHVFRKEQSRNM
KMERNHLGATLNECNQCFKVFSTKSSLTRHRKIHTGERPYGCSECGKSYSSRSYLAVHKRIHNGEKPYEC
NDCGKTFSSRSYLTVHKRIHNGEKPYECSDCGKTFSSSYLRPHLRHTGEKPYKCNQCFREFRTQSIFT
RHKRVHTGEGHYVCNQCAGKAFGRSSLSHYSHIHTGEYPYCHDCGRTRRRSNLTQHIRTHTGEKPYTC
NECGKSFTNSFSLTIHRRIHNGEKSYECSDCGKSFNVLSSVKKHMRHTGKKPYECNYCGKSFTSNSYLS
VHTRMHNRM

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_024341

ORF Size: 1290 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_024341.3](#)

RefSeq Size: 6002 bp

RefSeq ORF: 1293 bp

Locus ID: 79230

UniProt ID: [Q8N988](#)


Cytogenetics: 19p13.2

Domains: KRAB, zf-C2H2

Protein Families: Transcription Factors

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG216440