

Product datasheet for **RG216489**

Securin 2 (PTTG2) (NM_006607) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Tag: TurboGFP
Symbol: Securin 2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG216489 representing NM_006607
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCTACTCTGATCTACGTTGATAAGGAAATTGGAGAACCAGGCACCCGTGTGGCTGCCAAGGATGTGC
TGAAGCTGGAGTCTAGACCTTCAATCAAAGCATTAGATGGGATATCTCAAGTTTAAACACGACGTTTGG
CAAAACATACGATGCTCCATCAGCCTTACCTAAAGCTACCAGAAAGGCTTTGGGCACTGTCAACAGAGCT
ACAGAAAAGTCAGTAAAGACCAATGGACCCAGAAAACAAAACAGCCAAGCTTTTCTGCCAAAAAGATGA
CCGAGAAGACTGTAAACAAAAGTTCTGTTCTGCCTCAGATGACGCCTATCCAGAAATAGAAAAATT
CTTTCCTTCAATCTTCTAGACTTTGAGAGTTTGACCTGCCTGAAGAGCGCCAGATTGCACACCTCCCC
TTGAGTGGAGTGCCTCTCATGATCCTTGATGAGGAGGGAGAGCTTGAAAAGCTGTTTCAGCTGGGCCCCC
CTTCACCTGTGAAAATGCCCTCTCCACCATGGAATGCAATCTGTTTGCAGTCTCCTCAAGCATTCTGT
CGACCTGGATGTTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG216489 representing NM_006607
Red=Cloning site Green=Tags(s)

MATLIYVDKEIGEPGTRVAADVLKLESRPSIKALDGISQVLTRRFKTYDAPSALPKATRKALGTVNRA
TEKSVKTNQPRKQKQPSFSAKKMTEKTVKTKSSVPASDDAYPEIEKFFPNLLDFESFDLPEERQIAHLP
LSGVPLMILDEEGELEKLFQLGPPSPVKMPSPWECNLFVSVFKHSVDPGC

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: NM_006607

ORF Size: 575 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_006607.1](#), [NP_006598.1](#)

RefSeq Size:	576 bp
RefSeq ORF:	576 bp
Locus ID:	10744
UniProt ID:	Q9NZH5
Cytogenetics:	4p14
Protein Families:	Druggable Genome
Protein Pathways:	Cell cycle, Oocyte meiosis

Product images:



Circular map for RG216489