

Product datasheet for **RG216491**

SCRN3 (NM_024583) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	SCRN3
Synonyms:	SES3
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG216491 representing NM_024583
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGAACCTTTTCTGTGACACTTTCGTGGCATTACCTCCAGCAACAGTCGATAACAGGATTATTTTGG
GAAAAAATTCAGATAGACTCTATGATGAAGTACAAGAGGTGGTTTATTTTCTGCTGTAGTTCATGATAA
CCTGGGAGAACGCTTAAAGTGTACATATAGAAAATTGATCAAGTTCCTGAAACATATGCTGTTGTCCTG
AGTCGCCACGCGTGGTTGTGGGGGCGAGAAATGGGAGCCAATGAGCATGGAGTTTGCATTGGGAATGAAG
CTGTATGGGGAAGAGAAGAAAGTTTGTGATGAAGAAGCACTATTAGGAATGGACCTTGTGAGACTTGGCCT
TGAAAGAGCTGATACAGCTGAAAAAGCCCTCAATGTCATTGTTGACTTACTAGAAAAATATGGCCAGGGT
GGAAATTGCACAGAGGGTAGAATGGTATTTAGCTATCACAACAGTTTCTGATAGCTGATAGGAATGAAG
CCTGGATTCTGGAGACTGCAGGGAAGTACTGGGCAGCAGAAAAAGTACAAGAGGGAGTTCGTAATATTTT
TAATCAACTTTCCATAACAACCAAGATTGCCCGGGAACACCCAGACATGAGAACTATGCTAAGCGGAAA
GGTTGGTGGGATGGTAAAAAGGAGTTTGATTTTGTGTCAGCATATTCCTATCTTGACACAGCCAAGATGA
TGACTTCATCAGGCAGATACTGTGAGGGCTACAAGCTTCTAAATAAGCACAAGGAAATATAACTTTTGA
AACATGATGGAATTTCTCGAGATAAACCAAGTGGCATTAAATATGGAGGGAGAATTCCTGACCACTGCA
AGCATGGTTTCTATTTTACCTCAAGACTCCAGCCTTCTTGCACTTCTTTACAGGGACTCCTGATC
CTGAGAGATCTGTTTTAAGCCTTTTCAATTTGTGCCACATATTTCACAACTATTGGATACCAAGTTCACC
AACATTTGAACCTGAAGATCTAGTTAAAAAGAAATCACATTTTAAGCCTGACAGAAGACACCCACTCTAC
CAAAAACATCAACAGGCATTGGAAGTAGTAAATAAATGAGGAAAAAGCCAAATAATGTTGGACAACA
TGAGGAAACTGGAGAAAGAACTATTCAAGAGATGGAATCAATCCTTCAAAACAAGCATCTTGATGTGGA
GAAAATTGTTAATCTCTTCTCAGTGACAAAAGATGAAATTCAAATTTATCAGTCAAATTTATCAGTC
AAAGTTAGTTCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG216491 representing NM_024583
Red=Cloning site Green=Tags(s)

MEPFSCDTFVALPPATVDNRIIFGKNSDRLYDEVQEVVYFPAVVHDLGERLKCTYIEIDQVPETYAVVL
SRPAWLWGAEMGANEHGVCIGNEAVWGREEVCDEEALLGMDLVRLGLERADTAEKALNVIVDLLEKYGG
GNCTEGRMVFSYHNSFLIADRNEAWILETAGKYWAAEKVQEGVRNISNQLSITTKIAREHPDMRNYAKRK
GWWDGKKEFDFAAAYSILDTAKMNTSSGRYCEGYKLLNKHGKNITFETMMEILRDKPSGINMEGEFLT
TASMVSILPQDSSLPCIHFFTGTPDPERSVFKPFIHVPHISQLLDTSSPTFEEDLVKKKSHFKPDRRHPLY
QKHQQALEVVNNNEEKAKIMLDNMRKLEKELFREMESILQNKHLDEKIVNLFPPQCTKDEIQIYQSNLSV
KVSS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_024583

ORF Size: 1272 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_024583.5](#)

RefSeq Size: 2497 bp


RefSeq ORF: 1275 bp

Locus ID: 79634

UniProt ID: [Q0VDG4](#)

Cytogenetics: 2q31.1

Product images:



Circular map for RG216491