

Product datasheet for **RG216530**

GPR89C (NM_001097616) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GPR89C (NM_001097616) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	GPR89C
Synonyms:	GPHR; GPR89; GPR89C; SH120; UNQ192
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG216530 representing NM_001097616 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGTTTTAGGGATCTTATCCATAGAACAGCTCATCAGCCGGTTGGTGTGATTGGAGTGACTCTCATGG
CTCTTCTTTCTGGATTTGGTGTGTCAACTGCCATACACTTACATGTCTTACTTCCTCAGGAATGTGAC
TGACACGGATATTCTAGCCCTGGAACGGCGACTGCTGCAAACCATGGATATGATCATAAGCAAAAAGAAA
AGGATGGCAATGGCACGGAGAACAATGTTCCAGAAGGGGGAAGTGCATAACAAACCATCAGGTTTCTGGG
GAATGATAAAAAGTGTTACCACTTCAGCATCAGGAAGTAAAATCTTACTCTTATTCAACAGGAAGTGGA
TGCTTTGGAAGAATTAAGCAGGCAGCTTTTCTGGAAACAGCTGATCTATATGCTACCAAGGAGAGAATA
GAATACTCCAAAACCTTCAAGGGGAAATATTTTAAATTTCTTGTTACTTTTTCTCTATTTACTGTGTTT
GGAAAATTTTCATGGCTACCATCAATATTGTTTTGATCGAGTTGGGAAAACGGATCCTGTCAAGAGG
CATTGAGATCACTGTGAATTATCTGGGAATCCAATTTGATGTGAAGTTTGGTCCCAACACATTTCTTC
ATTCTTGTTGGAATAATCATCGTCACATCCATCAGAGGATTGCTGATCACTCTTACCAAGTTCTTTATG
CCATCTCTAGCAGTAAGTCCTCCAATGTCATTGTCCTGCTATTAGCACAGATAATGGGCATGTACTTTGT
CTCCTCTGTGCTGCTGATCCGAATGAGTATGCCTTTAGAATACCGCACCATAATCACTGAAGTCCTTGA
GAAGTGCAGTTCACTTCTATCACCGTTGGTTTGATGTGATCTTCTGGTCAGCGCTCTCTAGCATAC
TCTTCTCTATTTGGCTCACAAACAGGCACCAGAGAAGCAAATGGCACCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >RG216530 representing NM_001097616
 Red=Cloning site Green=Tags(s)

MFLGILSIEQLISRVGVIGVTLMLLLSGFGAVNCPYTYMSYFLRNVTDDILALERRLLQTMDMIISKKK
 RMAMARRTMFQKGEVHNKPSGFWGMIKSVTTSASGSENLTLIQQEVDALEELSRQLFLETADLYATKERI
 EYSKTFKGKYFNFLGYFFSIYCVWKIFMATINIVFDRVGKTDVPTRGIEITVNYLGIQFDVKFWSQHISF
 ILVGIIIVTSIRGLLITLTKFFYAISSSKSSNVIVLLLAQIMGMVFVSSVLLIRMSMPLEYRTIITEVLG
 ELQFNFYHRWFDVIFLVLSALSSILFLYLAHKQAPEKQMAP

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001097616

ORF Size: 960 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001097616.1, NP_001091085.1

RefSeq Size: 1569 bp


RefSeq ORF: 962 bp

Locus ID: 728932

Cytogenetics: 1q21.1

Protein Families: Transmembrane

Product images:



Circular map for RG216530