

Product datasheet for **RG216543**

C4orf33 (NM_001099783) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C4orf33 (NM_001099783) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: C4orf33
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG216543 representing NM_001099783
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGATTTTAAAATTGAACACACTTGGGATGGTTTTCCAGTGAAGCATGAGCCCGTATTTATCAGGCTGA
 ATCCAGGTGACAGAGGAGTGATGATGGACATTAGTGCTCCATTTTTTCAGGGATCCTCCAGCCCCACTTGG
 AGAACCAGGAAAACCTTTCAATGAAGTGTGGGATTATGAAGTTGTGGAAGCATTTTTCTTGAATGATATA
 ACTGAGCAATATTTAGAAGTTGAAGTGTCCACGGACAGCATTAGTGCTTTTACTTTCTGGAAGAA
 GAAATGTGTGGAACAAGAAGTTCCTTTATCGTTCAGAAATGTCCAGAGGAGAGACAAAATGGGAAGGCAA
 AGCTTATCTCCCTTGGAGTTATTTTCCACCAAATGTGACAAAATTCATTTCAATTTGCAATTCATGGATCA
 AAAGATAAACGAAGTTATGAAGCTCTTTACCCTGTACCTCAGCATGAACTGCAGCAAGGACAAAAACCTG
 ATTTCCATTGCCTAGAATACTTCAAGTCTTTCAATTTTAAACACTGCTTGGAGAAGAGTGGAAACAAC
 GGAATCAGACCTGTGGCTAATAGAGAAATGTGATATA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG216543 representing NM_001099783
 Red=Cloning site Green=Tags(s)

MDFKIEHTWDGFPVKHEPVFIRLNPGRGVMMDISAPFFRDPPAPLGEKGFNFELWDYEVVEAFFLNDI
 TEQYLEVELCPHGQHLVLLL SGRNVWKQELPLSFRMSRGETKWEKAYLPWSYFPPNVTKFNSFAIHGS
 KDKRSYEALYPVQHELQQGQKPDFHCLEYFKSFNFNTLLGEEWKQPESDLWLIKCDI

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001099783

ORF Size: 597 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001099783.2</u>
RefSeq Size:	1601 bp
RefSeq ORF:	600 bp
Locus ID:	132321
UniProt ID:	<u>Q8N1A6</u>
Cytogenetics:	4q28.2