

Product datasheet for **RG216576**

ZNF641 (NM_152320) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	ZNF641
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG216576 representing NM_152320
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCAAGCTGAGGACAGGTCACAATTTGGGTCTGCAGCAGAAATGCTTTCAGAACAGACAGCAGCGCTGG
GGACAGGATGGGAATCAATGAATGTACAGCTGGATGGAGCAGAGCCCCAGGTGGAAAGGGGAAGCCAGGA
AGAGCGGCCATGGAGAACAGTACCAGGACCTCTGGAGCACCTGTGCTGTGACCTTGAAGAGGAGCCACAG
TCCCTTCAGGAGAAGGCTCAGTCTGCTCCCTGGGTTCTCTGCAATCCCCAGGAGGGGAACACTGGAGACT
GGGAGATGGCAGCTGCACTTCTTGGGCTGGATCACAGGCGCTGGTAACCATCAAGGATGTGTCACTGTG
CTTCTCTCAGGAGGAGTGGCGGAGCCTGGACCCCTCTCAGACAGACTTTTATGGAGAATATGTCATGCAG
GAAACTGTGGGATAGTAGTCTCTCTGAGATTTCCAATTCCTCAAACTGGACATGCTTTCTCAACTAGAAG
GAGGAGAAGAACAATGGGTCCCTGACCCCGAGGACTTAGAGGAGAGGGACATTCTGAGGGTCACATATAC
AGGAGGTGGAAGTGAACATGAGGGGGATACCCCTGAACTAGAAGCAGAACCTCCAGAATGTTATCCAGC
GTGTCTGAAGATACTGTTCTCTGGAACCCGGAGCATGATGAGAGCTGGGATTCCATGCCAGCAGCTCCA
GAGGAATGCTCCTGGGGCCCCCTTTCTTCAGGAAGATAGTTTCTCAAACCTGCTGTGTAGCACAGAGAT
GGATTCCCTGTTAAGACCCACACATGCCCCAGTGTGGGAAACAGTTTGTATGGGGTTCACCTTGCC
AGGCATCAACAAACACACACTGGGGAGAGACCTACAGCTGCCTCAAGTGTGAGAAGACCTTTGGGCGAA
GACATCACCTCATCAGGCACCAGAAAACCCACCTACATGACAAGACCAGCAGGTGCTCTGAGTGTGGTAA
GAATTTCCGATGCAACTCCCATCTGGCCAGCCACCAGAGAGTGCATGCAGAAGGCAAACTCTGCAAAAGGC
CAAGAGGTTGGAGAGAGCCCTGGCACAAGGAAACGGCAGCGTGGCCCAAGTGCCTGCAAAAGTGTACAGTGT
GCACCTGAATGTGGGAAGAGCTTTGGCCGAAGGCACCACCTTGTGAGACACTGGTGACCCACACTGGGGA
GAAGCCCTTCAGTGCCCTCGCTGTGAGAAGAGCTTTGGCCGAAACATCACCTGGACAGGCACCTGCTC
ACCCACCAGGGACAAAGTCCCGGAACAGCTGGGACAGAGGAACATCTGTCTTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence:

>RG216576 representing NM_152320

Red=Cloning site Green=Tags(s)


MQAEDRSQFGSAAEMLSEQTAALGTGWESMNVQLDGAEPQVERGSQEERPWRTVPGPLEHLCCDLEEEPQ
SLQEKAQSAPWVPAIPQEGNTGDWEMAAALLAAGSQGLVTIKDVSLCFSQEEWRSLDPSQTDYFGEYVMQ
ENCGIIVSLRFPPIPKLDMLSQLEGGEQWVPDPQDLEERDILRVITYTGGGSEHEGDTPELEAPPRMLSS
VSEDTVLWNPEHDESWDSMPSSSRGMLLGPPFLQEDSFSNLLCSTEMDSLRLPHTCPQCGKQFVWGSHLA
RHQQTHTGERPYSLCKCEKTFGRRHHLIRHQKTHLHDKTSRCSECKNFRCNHSLASHQRVHAEGKSCKG
QEVGESPGTRKRQRAPPVPKCHVCTECGKSFGRRHHLVRHWLTHTEKPFQCPRCEKSFGRKHHLDRHLL
THQGQSPRNSWDRGTSVF

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_152320

ORF Size:

1314 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_152320.1](#), [NP_689533.1](#)

RefSeq Size: 2326 bp

RefSeq ORF: 1317 bp

Locus ID: 121274

UniProt ID: [Q96N77](#)


Cytogenetics: 12q13.11

Domains: KRAB, zf-C2H2

Protein Families: Transcription Factors

Gene Summary: Transcriptional activator. Activates transcriptional activities of SRE and AP-1.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG216576