

## Product datasheet for **RG216609**

### DSN1 (NM\_024918) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	DSN1 (NM_024918) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	DSN1
Synonyms:	C20orf172; dj469A13.2; hKNL-3; KNL3; MIS13
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG216609 representing NM_024918 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGACTTCAGTGACTAGATCAGAGATCATAGATGAAAAAGGACCAGTGATGTCTAAGACTCATGATCATC  
AATTGGAATCAAGTCTCAGTCCTGTGGAAGTGTTGCTAAAACATCTGCCTCCCTGGAGATGAATCAAGG  
CGTTTCAGAGGAAAGAATTCACCTTGGCTCTAGCCCTAAAAAGGGGAAATTGTGATCTCAGCCACCAG  
GAAAGACTTCAGTCGAAGTCCCTTCATTTGTCTCCTCAAGAACAATCTGCCAGTTATCAAGACAGGAGGC  
AATCCTGGCGGCGAGCAAGTATGAAAGAAACGAACCGGCGGAAGTCGCTGCATCCCATCCAGGGCAT  
CACAGAGCTCAGCCGGTCTATCAGTGTCGATTTAGCAGAAAGCAAACGGCTTGGCTGTCTCCTGCTTCC  
AGTTCCAGTTCTCTATTAGAACTTGAACCTTCTAAGGGACACTAAGGGCTTCAGTCTTGAAGTT  
TTAGAGCCAAAGCATCTTCTCTTCTGAAGAATTGAAACATTTTGCAGACGGACTGGAACTGATGGAAC  
TCTACAAAAATGTTTTGAAGATTCAAATGGAAAAGCATCAGATTTTTCTTGGAAAGCATCTGTGGCTGAG  
ATGAAGGAATACATAACAAAGTTTTCTTTAGAACGTCAGACTTGGGATCAGCTCTTGCTTCACTACCAGC  
AGGAGGCTAAAGAGATATTGTCAGAGGATCAACTGAGGCCAAAATTACTGAGGTCAAAGTGAACCTAT  
GACATATCTTGGTCTTCTCAGAATGAAGTTCTTAATACAAAACCTGACTACCAGAAAATATTACAGAAC  
CAGAGCAAAGTCTTTGACTGTATGGAGTTGGTGATGGAAGTCAAGGATCAGTGAACAGCTGCAGG  
CCTTTATGGATGAAAGTACCCAGTGCTTCCAGAAGGTGTGATCAGCTCGGAAAAGAGAAGCATGCAACA  
ATTAGATCCCTCACCAGCTCGAAAACCTGTTGAAGCTTCAGCTACAGAACCACCTGCCATACATGGATCT  
GGATCTGGATCTTGTGAG

**ACGGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG216609 representing NM\_024918  
 Red=Cloning site Green=Tags(s)

MTSVTRSEIIDEKGPVMSKTHDHQLESSLSPVEVFAKTSASLEMNQGVSEERIHGSSPKKGGNCDLSHQ  
 ERLQSKSLHLSPQEQSASYQDRRQSWRRASMKETNRRKSLHPIHQGITELSRISVDLAESKRLGCLLLS  
 SSQFSIQKLEPFLRDTKGFSLSEFRKASSLSEELKHFADGLETDLQKCFEDSNGKASDFLEASVAE  
 MKEYITKFSLERQTWDQLLLHYQQEAKEILSRGSTEAKITEVKVEPMTYLGSSQNEVLNTPDYQKILQN  
 QSKVFDCEMELVMDELQGSVKQLQAFMDESTQCFQKVSVQLGKRSMQQLDPSARKLLKLQLQNPPAIHGS  
 GSGSCQ

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_024918

**ORF Size:** 1068 bp


**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
- RefSeq:** [NM\\_024918.2](#), [NP\\_079194.2](#)
- RefSeq Size:** 2095 bp
- RefSeq ORF:** 1071 bp
- Locus ID:** 79980
- UniProt ID:** [Q9H410](#)
- Cytogenetics:** 20q11.23
- Gene Summary:** This gene encodes a kinetochore protein that functions as part of the minichromosome instability-12 centromere complex. The encoded protein is required for proper kinetochore assembly and progression through the cell cycle. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2009]

**Product images:**



Circular map for RG216609