

Product datasheet for **RG216631**

C14orf184 (NM_001080113) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Tag: TurboGFP
Symbol: C14orf184
Synonyms: FLJ33842; PP8961
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG216631 representing NM_001080113
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGAGTGTCCCTAAATTAGTCAACCATGGCTGGGAGACAGAGCCATGCAATCCAGACCCAGCTCTTCTGC
TCACTGGGGCCGTTCCCGGAGAAGGAGAAACACGCTGGAGCGGGCACTCCACGAGTTGTCTGGGGAACA
AGAACTCGGATCCTGGCGGGAGAGAGCGAACAGTACACCCGAGCGCGGCTGTACCCGCTTTGAAACAC
GTAAGTGCCCGCCAACCAGCATCCAGAAACGCCTGGGAGCGGGCGCGCGGACGAGGAGAGACCGCG
CTGGGAGGGCGCGGAAGGCGAGGGCGCTGCCGCTCCACGCTGAGCTCACTGCGCGGGGCGGCCCGCGA
GGCAGGGCGCGGGCGGAGGCGCCCGAGTGGTCCGTCTCGGAGGCGCAGACGCGGTTCTCTCGAGGGC
GCTCCCGCTGTGAGGAGCAAATCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG216631 representing NM_001080113
Red=Cloning site Green=Tags(s)


MSVPKLVNHGWETPCNPDPALLLTGAVPGEGETTLERLPRVVGEQELGSWRERANSTPERRLSPALKH
VSARPTSIPETPGERGARTRRDRAGRARKGEGAAASTLSSLRGAAPAEAGRAEAPWSVSGGADAVPLEG
APAVRSKS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: NM_001080113

ORF Size: 444 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_001080113.2](#), [NP_001073582.2](#)

RefSeq Size: 2596 bp
RefSeq ORF: 446 bp
Locus ID: 650662

Cytogenetics: 14q32.12

Gene Summary: DISCONTINUED: This record has been withdrawn by NCBI staff. This record was defined by AF318368.1 and AK091161.1, which appear to be cloning artifacts and are not sufficient evidence to define a distinct locus.

Product images:



Circular map for RG216631