

Product datasheet for **RG216665**

MASP2 (NM_006610) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MASP2 (NM_006610) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MASP2
Synonyms:	MAP-2; MAP19; MASP-2; MASP1P1; sMAP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RG216665 representing NM_006610
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAGGCTGCTGACCCTCTGGCCTTCTGTGTGGCTCGGTGGCCACCCCTTGGGCCCGAAGTGGCCTG
 AACCTGTGTTTCGGGCGCCTGGCATCCCCGGCTTCCAGGGGAGTATGCCAATGACCAGGAGCGGCCTG
 GACCCTGACTGCACCCCCGGCTACCGCTGCGCCTCTACTTCACCCACTTCGACCTGGAGCTCTCCAC
 CTCTGCGAGTACGACTTCGTCAAGCTGAGCTCGGGGGCCAAGGTCTGGCCACGCTGTGCGGGCAGGAGA
 GCACAGACACGGAGCGGGCCCTGGCAAGGACACTTTCTACTCGTGGGCTCCAGCTGGACATTACCTT
 CCGCTCCGACTACTCCAACGAGAAGCCGTTACGGGGTTCGAGGCCTTCTATGCAGCCGAGGACATTGAC
 GAGTGCCAGGTGGCCCCGGGAGAGGCGCCACCTGCGACCACCTGCCACAACCACCTGGGCGGTTTCT
 ACTGCTCTGCCGCGCAGGCTACGTCTGCACCCTAACAAAGCGCACCTGCTCAGCCCTGTCTCCGGCCA
 GGTCTTACCCAGAGGTCTGGGGAGCTCAGCAGCCCTGAATACCCACGGCCGTATCCCAAATCTCCAGT
 TGCACTTACAGCATCAGCCTGGAGGAGGGTTCAGTGTCACTTCTGGACTTTGTGGAGTCTTCGATGTGG
 AGACACACCCTGAAACCCTGTGTCCCTACGACTTCTCAAGATCAAACAGACAGAGAAGAACATGGCCC
 ATTCTGTGGGAAGACATTGCCCCACAGGATTGAAACAAAAGCAACACGGTGACCATCACCTTTGTCA
 GATGAATCAGGAGACCACACAGGCTGGAAGATCCACTACAGCAGCAGCGCAGCCTTGCCCTTATCCGA
 TGGCGCCACCTAATGGCCACGTTTACCTGTGCAAGCCAAATACATCCTGAAAGACAGCTTCTCCATCTT
 TTGCGAGACTGGCTATGAGCTTCTGCAAGGTCACTTGGCCCTGAAATCCTTTACTGCAGTTTGTGAGAAA
 GATGGATCTTGGGACCGCCAATGCCCGCTGCAGCATTGTTGACTGTGGCCCTCCTGATGATCTACCCA
 AGAGACCTTCTACACAATGAAAGTGAATGATGGTAAATATGTGTGTGAGGCTGATGGATTCTGGACGAGC
 TCCAAAGGAGAAAAATCACTCCAGTCTGTGAGCCTGTTTGTGGACTATCAGCCCGCACAACAGGAGGGC
 GTATATATGGAGGGCAAAGGCCAAAACCTGGTGATTTTCTTGGCAAGTCTGATATTAGGTGGAACCAC
 AGCAGCAGGTGCACTTTTATATGACAACTGGGTCTAACAGCTGCTCATGCCGTCTATGAGCAAAAACAT
 GATGCATCCGCCCTGGACATTGGAATGGGCACCCTGAAAAGACTATCACCTCATTATACACAAGCCTGGT
 CTGAAGCTGTTTTATACATGAAGTTTACTCATGATGCTGGCTTTGACAATGACATAGCACTGATTAA
 ATTGAATAACAAAGTTGTAATCAATAGCAACATCACGCCATTTTGTCTGCCAAGAAAAGAAGCTGAATCC
 TTTATGAGGACAGATGACATTGGAAGTGCATCTGGATGGGGATTAAACCAAAGGGGTTTTCTTGCTAGAA
 ATCTAATGTATGTCGACATACCGATTGTTGACCATCAAAAATGACTGCTGCATATGAAAAGCCACCCTA
 TCCAAGGGGAAGTGAACCTGCTAACATGCTTTGTGCTGGCTTAGAAAAGTGGGGGCAAGGACAGCTGCAGA
 GGTGACAGCGGAGGGGCACTGGTGTCTAGATAGTGAACAGAGAGGTGGTTTGTGGGAGGAATAGTGT
 CCTGGGTTCCATGAATTGTGGGAAGCAGGTGAGTACACAAAAGTTATTAACATATATCC
 CTGGATCGAGAACATAATTAGTGATTTT

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA


Protein Sequence: >RG216665 representing NM_006610
 Red=Cloning site Green=Tags(s)

```
MRLLTLLGLLCGSVATPLGPKWPEPVFGRLASPGFPGEYANDQERRWLTAPPGYRLRLYFTHFDLELSH
LCEYDFVKLSSGAKVLATLCGQESTDTERAPGKDTFYSLGSSLDITFRSDYSNEKPF TGFEAFYAAEDID
ECQVAPGEAPTCDHHCHNHLGGFYCSCRAGYVLRNKRKTCALCSGQVFTQRSGELSSPEYPRYPKLLS
CTYSISLEEGFSVILDFVESFDVETHPETLCPYDFLKIQTDRREEHGPF CGKTLPHRIETKSNTVTITFVT
DESGDHTGWKIHYTSTAQPCPYMAPPNGHVSPVQAKYILKDSFSIFCETGYELLQGHLPKLSFTAVCQK
DGSWDRMPACSI VDCGPPDDLPSGRVEYITGPGVTTYKAVIQYSCEETFYTMKVNDGKYVCEADGFWTS
SKGEKSLPVCEPVCGLSARTTGRIYGGQKAKPGDFPWQVILGGTTAAGALLYDNWVLTAAHAVYEQKH
DASALDIRMGT LKRLSPHYTQAWSEAVF IHEGYTHDAGFDNDIALIKLNKVVINSNITPICLPRKEAES
FMRTDDIGTAGSGWGLTQRGFLARNLMYVDIPIVDHQKCTAAYEKPPYPRGSVTANMLCAGLESGGKDCSR
GDSGGALVFLDSETERWFVGGIVSWGSMNCGEAGQYGVYTKVINIPIWIENIISDF
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_006610

ORF Size: 2058 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_006610.4
RefSeq Size:	2460 bp
RefSeq ORF:	2061 bp
Locus ID:	10747
UniProt ID:	O00187
Cytogenetics:	1p36.22
Domains:	CCP, CUB, Tryp_SPc, EGF_CA, EGF
Protein Families:	Druggable Genome, Protease, Secreted Protein
Protein Pathways:	Complement and coagulation cascades
Gene Summary:	<p>This gene encodes a member of the peptidase S1 family of serine proteases. The encoded preproprotein is proteolytically processed to generate A and B chains that heterodimerize to form the mature protease. This protease cleaves complement components C2 and C4 in order to generate C3 convertase in the lectin pathway of the complement system. The encoded protease also plays a role in the coagulation cascade through cleavage of prothrombin to form thrombin. Myocardial infarction and acute stroke patients exhibit reduced serum concentrations of the encoded protein. Alternative splicing results in multiple transcript variants, at least one of which encodes an isoform that is proteolytically processed. [provided by RefSeq, Feb 2016]</p>

Product images:



Circular map for RG216665