

Product datasheet for **RG216764**

PRCD (NM_001077620) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Tag: TurboGFP

Symbol: PRCD

Synonyms: RP36

Mammalian Cell Selection: Neomycin

Vector: pCMV6-AC-GFP (PS100010)

E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG216764 representing NM_001077620
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTGCACCACCCTTTTCCTGCTCAGCACCTGGCCATGCTCTGGCGCCGCCGATTTGCCAACCGAGTCC
AACCAGAGCCGAGCGACGTGGATGGGGCAGCTAGGGGCAGCAGCTTGGATGCGGACCCTCAGTCCTCAGG
CAGGGAGAAAGAACCTCTGAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG216764 representing NM_001077620
Red=Cloning site Green=Tags(s)

MCTTLFLLSTLAMLWRRRFANRVQPEPSDVGAAARGSSLDADPQSSGREKEPLK


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



Cloning Scheme:

Cloning sites used for ORF Shutting:



EcoRI BamHI KpnI RBS Kozac Consensus SgfI AscI
 CTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGATCGCCGCCGCGCAGATCT
 HindIII NheI RsrII MluI NotI XhoI GFP Tag
 CAAGCTTAAGTCTAGTAGCGGACCG ACG CGT ACG CGG CCG CTC GAG ATG GAG AGC GAC --- --- ---
 T R T R P L E M E S D - - -
 --- --- GAA GAA AGA GTT TAA ACGGCCGCCGCGGAGCT
 - - E E R V Stop

ACCN: NM_001077620

ORF Size: 162 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001077620.3</u>
RefSeq Size:	937 bp
RefSeq ORF:	165 bp
Locus ID:	768206
UniProt ID:	<u>Q00LT1</u>
Cytogenetics:	17q25.1
Gene Summary:	This gene is predominantly expressed in the retina, and mutations in this gene are the cause of autosomal recessive retinal degeneration in both humans and dogs. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Mar 2010]

Product images:



Circular map for RG216764