

## Product datasheet for **RG216769**

### **C9orf5 (TMEM245) (NM\_032012) Human Tagged ORF Clone**

#### **Product data:**

|                                  |                         |
|----------------------------------|-------------------------|
| <b>Product Type:</b>             | Expression Plasmids     |
| <b>Tag:</b>                      | TurboGFP                |
| <b>Symbol:</b>                   | C9orf5                  |
| <b>Synonyms:</b>                 | C9orf5; CG-2; CG2       |
| <b>Mammalian Cell Selection:</b> | Neomycin                |
| <b>Vector:</b>                   | pCMV6-AC-GFP (PS100010) |
| <b>E. coli Selection:</b>        | Ampicillin (100 ug/mL)  |



[View online »](#)

ORF Nucleotide Sequence: >RG216769 representing NM\_032012  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCGCGATCGCC

ATGCCGACGGCGGCCCTAAGGACGCGCAAGCCTGCGGAGCTCTCCGGGCGGCCGCGGGTCC  
 CGCGCGCGGTTCGGGCCGAGTGGCGGTGGCGGGAGACCCCGGACCGCGCGCTGGCGCTGCGCTTCGA  
 CAAGCCCATTAAGCAGGCTTCTACAACACCGGGCCGTGCTGTTCTGCTGCTGCTGCTGCGGCGCGCG  
 GTGCTGGTCTACTTCACTCTGGAGGCTTCTGCGGCCGCTGCTCTGGGCCGTGCTATGCGGCACTTTTC  
 TGCACCCCTTCAAGAGCTCGCTGACGCGCTGGGCCGCACTGGCTGCAGCGCTGCACCGCGCGCACAC  
 GCCCATCGTCTGGCCGCGCTGCTCTGCGCTGCTGCTTCGTCGACTACGGCGTCGAGGCCCTGGGCGAG  
 CAGGCGCTGCGCCGCCGCGCTGCTCTGCTGCTGCGCGCGGCCGCGCTCTGTACGGCTCTACT  
 GCCTCGGCAGCTACTTGGCGTGCAGGTGCTGCTGGTGCACGCTGCCACGCTCATCTGCCGCGGGTGA  
 CTACTTCAGCAGCTGTGGATCTGGACGTTGGTGGTGGCTATGTGTTGACTGTTTCATTCAAGTGAAT  
 GCAAGCACTGAACGCTACTTGAGAGCAGTTCAATTCTGTCTGGATTATATTGCTTTTTCATTTAGCAT  
 CCCTGGCTGGCTCATGGAGAATCCAGTATTCTGGTTATTGTTTTCTGATGTCTGTGGGTACCTCTA  
 TGAAAAACAGAAATGAAAAGAGTCTTCTGGGCGAGATTACAGGGCAGGTTATCTCCATGGCAGCTTCT  
 ACTCTGGCAAATTTGGCCATCTCTATCACAGGATGAATCAAGCAGTGAAGACCAGCCCTCCACTCAGC  
 CTGCAGAAGCAGTGGACAGGGGAGAATCCGCTCCAACGTTGTCCACCTCCCTTCACCTCCTCCCTTC  
 ACCCACTTCCCTTCACCTACTCTGGGCGAGCAAGGCTGAAATAGGAACGTTTCTAGAAAGAAGAAA  
 ACTAGTGACATCTACTTTGTGTCTAGTTTGGCCATTGTCGTCATGCAGATCTGGTTGAACCTGTGGA  
 TTGTGCAGTTGCTGCCAGTCCGATTGCAGTCTGGATACTCAAAAGCTTGTCACTTCTTTGGAGTTGT  
 GGATTTCTAGAGAAACGCTACCATGTGTGGTGGGCGATTATAGAAAGCTTCTGAAGGAGCGCGAGGGA  
 GCTCTCGCGCTTGGCCATTGTCGGGCTTGAAAAATCTTGTTAAAGTTGATAGCAAGCTCTGGCACT  
 GGCTAAATAAGAAGATGATCATCTGGTTGGAGAAGATGTAGATAAAATAATTAGCATTTTCATCATATT  
 TTTGTTAGTGATAGGAATCTTCTTTAGCCCTACTCCTGACTGCAAGGTACATCAAGAGAGTGATACAC  
 ATGATTGAAGTCACAAGTAATTTGATTAAAGAACTCTAGCAAAATCACCTGAGTGGGCAAAATGGCTTC  
 CTGAGGCTCAGGTAGTCCAAAGAGCCCTGAATTCTGCGGCTAACAACGTGTATCAGTATGGACGAGAATG  
 GATAACTACAAGCTCCATAAAATTTAGGAGATAAGGTGAACAATACTGCTGTAATTGAAAAGCAAGTA  
 CTAGAACTTTGGGACAGACTGTACTCTTGGTTTGTAAAGAATGAACACACTCTGGAAGGCACAAAG  
 GACAGAAGTTGCATGTCAGTCGTCAGAAATAGCTGGCTGGGAGACATTCTGGACTGGCAGGATATTGTTTC  
 CTTTGTTCACGAGAACATTGAGACATTTCTTCGATCTTGGAGTCTCTGTGGATCGTTATGAGCCGGAAT  
 GTGAGCCTGCTGTTACCACCTGTCACTACACTCTTGACCATCCTCTTACAGCGGGACAGCCCTTCTCA  
 ATTTTGTACTCTCTGATAATTTTCTGACCACACTATTTATCTATTAAGTTCCAGTGATGAGTACTA  
 CAAGCCAGTGAAGTGGGTGATAAGCCTGACTCCACTATCTCAGCCAGGTCCTTCTTCTAATATTATTGGC  
 CAGTCTGTGGAAGAAGCTATCAGAGGGGTGTTGATGCTTCCCTCAAAATGGCTGGCTTCTATGGATTGT  
 ATACCTGGCTGACTCATACTATGTTTGGCATCAATATTGTCTTCAACCATCAGCATTAGCAGCAATCCT  
 TGGAGCAGTGCCATTCTGGGACATACTGGGACGAGTACCTGCAGTTCTTGACCTGTGGCTGACACAA  
 GGGTTAGGATGCAAGGCCATTTTACTGTTGATTTTTCATCTCTTGCCAACATACTTTGTAGATACTGCAA  
 TCTACTCTGACATATCAGGAGGTGGCCATCCTTACCTGACAGGCTTGGCAGTGGCCGGTGGAGCATACTA  
 CCTAGGCCCTGGAAGGAGCAATCATCGGTCTATTCTTCTCTGCATACTTGTGTTGCTTCCAATATCTAT  
 AGTGCCATGCTAGTGAGTCCCACGAATTCAGTTCACGCCAAACCAGACCCCATGGCTGCTCAGCCTC  
 AGCGGACTTTCGCTGACATTTCTGAAGATCTGAAATCTTCAGTAGGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:**

>RG216769 representing NM\_032012

Red=Cloning site Green=Tags(s)


```
MADGGGPKDAPSLRSSPGPAPRVPRVAVGPGGGGETPRTAALALRFDKPIKQAFYNTGAVLFVCLCCGAA
VLVYFILEAFLRPLLWAVLCGTFLLHPFKSSLTRLGRHWLQRLHRAHTPIVLAALLPLCFVDYGVGEALGE
QALRRRRLLLLLGAGGPLLYGLYCLGSYLGQVLLVHAATLICRGLDYFSSLWIWTLVVGYYLTVSFKWN
ASTERYLRAVSIPVWIILLFHLASLAGSWRIPVFLVIVFLMSVGTLYEKQNGKESSGAELPGQVISMAAS
TLANLAISITGYESSSEDQPSTQPAEAVDRGESAPTLSTSPSPSSPSPTSPSPTLGRRRPEIGTFLRKKK
TSDIYFVSLVWAIIVMQIWLNLWIVQLLPVPIAVWILKKLVHFGVVDFFLEKRYHVWVGIIIESFLKERQG
ALAPWPIVGLGKFLKVDKSLWHWLNKKMIWLEKMLDKIISIFIIIFLLVIGTLLALLLTAKVHQESVH
MIEVTSNLINETLANHPEWANWLPEAQVVQALNSAANNVYQYGREWITKHLKILGDKVNNTAVIEKQV
LELWDRLYHSWFVKNVTHSGRHKGQKLHVSQRNSWLGDILDWQDIVSFVHENIETFLSILESLWIVMSRN
VSLLFTTVTLLTILFYSGTALLNFVLSLIIFLTTLFYLLSSSDEYYKPVKWISLTPLSQPGPSSNIIG
QSVEEAIRGVFDASLKMGAFYGLYTWLTHTMFGINIVFIPSALAAILGAVPFLGTYYAAVPAVLDLWLTQ
GLGCKAIIILLIFHLLPTYFVDTAIYSDISGGGHPYLTGLAVAGGAYYLGLEGAIGPILLCILVVASNIY
SAMLVSPTNSVPTPNQTPWPAQPQRTFRDISEDLKSSVG
```

TRTRPLE - GFP Tag - V

**Restriction Sites:**

Sgfl-MluI

## Cloning Scheme:



**ACCN:** NM\_032012

**ORF Size:** 2637 bp


**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>Note:</b>                  | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>RefSeq:</b>                | <u>NM_032012.4</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>RefSeq Size:</b>           | 8006 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 2640 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>              | 23731                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>UniProt ID:</b>            | <u>Q9H330</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Cytogenetics:</b>          | 9q31.3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Protein Families:</b>      | Transmembrane                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |

**Product images:**



Circular map for RG216769